

# Package ‘scpdata’

April 30, 2026

**Type** Package

**Title** Single-Cell Proteomics Data Package

**Version** 1.21.0

**Description** The package disseminates mass spectrometry (MS)-based single-cell proteomics (SCP) datasets. The data were collected from published work and formatted using the `scp` data structure. The data sets contain quantitative information at spectrum, peptide and/or protein level for single cells or minute sample amounts.

**Depends** R (>= 4.2.0), QFeatures, ExperimentHub

**Imports** utils, AnnotationHub, SingleCellExperiment, S4Vectors

**Suggests** scp, magrittr, dplyr, knitr, BiocStyle, BiocCheck, rmarkdown, testthat

**biocViews** ExperimentData, ExpressionData, ExperimentHub, ReproducibleResearch, MassSpectrometryData, Proteome, SingleCellData, PackageTypeData

**License** GPL-2

**Encoding** UTF-8

**LazyData** false

**VignetteBuilder** knitr

**Roxygen** list(markdown = TRUE)

**RoxygenNote** 7.3.3

**git\_url** <https://git.bioconductor.org/packages/scpdata>

**git\_branch** devel

**git\_last\_commit** dec730b

**git\_last\_commit\_date** 2026-04-28

**Repository** Bioconductor 3.24

**Date/Publication** 2026-04-30

**Author** Christophe Vanderaa [aut, cre] (ORCID: <https://orcid.org/0000-0001-7443-5427>),  
Laurent Gatto [aut] (ORCID: <https://orcid.org/0000-0002-1520-2268>),  
Enes Ayar [dtc] (ORCID: <https://orcid.org/0000-0001-8881-9800>),  
Samuel Grégoire [dtc] (ORCID: <https://orcid.org/0000-0002-0900-7677>)

**Maintainer** Christophe Vanderaa <[vanderaa.christophe@gmail.com](mailto:vanderaa.christophe@gmail.com)>

## Contents

ai2025a	2
brunner2022	3
cong2020AC	5
derks2022	6
dou2019_boosting	8
dou2019_lysatess	10
dou2019_mouse	12
gregoire2023_mixCTRL	14
guise2024	16
hu2023_K562	18
khan2023	20
krull2024	22
leduc2022	23
leduc2022_plexDIA	24
leduc2022_pSCoPE	26
liang2020_hela	28
petrosius2023_AstralAML	29
petrosius2023_mES	31
schoof2021	33
scpdata	35
specht2019v2	36
specht2019v3	38
williams2020_lfq	40
williams2020_tmt	42
woo2022_lung	43
woo2022_macrophage	45
wu2026	46
zhu2018MCP	48
zhu2018NC_hela	49
zhu2018NC_islets	51
zhu2018NC_lysatess	52
zhu2019EL	54
<b>Index</b>	<b>56</b>

---

 ai2025a

*Ai et al., Cardiomyocytes Single Cell Proteomics, MCP (2025)*


---

### Description

The manuscript analysed cardiomyocytes derived from iPSC cells and adult hearts. This annotated QFeatures objects contains the complete data, as analysed with the scp package. The re-analysis was performed as part of the EuBIC 2025 developer meeting hackathon dedicated to single-cell proteomics data (<https://lgiatto.github.io/2025-EuBIC-SCP-hackathon/>).

The data contains a total of 304 assay. The 299 first assays contain the individual single cardiomyocytes acquisitions. The following assays contain the joined precursor data, the log-transformed precursor data, the peptide data, the protein data, and the scplainer-modelled precursor data (see the link above for details).

The `QFeatures` object was generated from the DIA-NN report file, available from the `MsDataHub` package (see `?MsDataHub:ai2025_aCMS_report.tsv`). The report file is downloaded from the MassIVE MSV000094438 (doi:10.25345/C5T727S7Q) repository.

### Usage

```
ai2025a
```

### Format

An object of class function of length 1.

### References

- Ai, Lizhuo, Aleksandra Binek, Vladimir Zhemkov, Jae Hyung Cho, Ali Haghani, Simion Kreimer, Edo Israely, et al. 2025. “Single Cell Proteomics Reveals Specific Cellular Subtypes in Cardiomyocytes Derived from Human iPSCs and Adult Hearts.” *Molecular & Cellular Proteomics: MCP*, no. 100910 (January): 100910.

### Examples

```
ai2025a()
```

---

brunner2022

*Brunner et al. 2022 (Mol. Syst. Biol.): cell cycle state study*

---

### Description

Single cell proteomics data acquired by the Mann Lab using a newly designed timsTOF instrument, referred to as timsTOF-SCP. The dataset contains quantitative information from single-cells blocked at 4 cell cycle stages: G1, G1-S, G2, G2-M. The data was acquired using a label-free sample preparation protocol combined to a data independent (DIA) acquisition mode.

### Usage

```
brunner2022
```

### Format

A `QFeatures::QFeatures` object with 435 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object.

- Assay 1-434: DIA-NN main output report table split for each acquisition run. Since each run acquires 1 single cell, each assay contains a single column. It contains the results of the spectrum identification and quantification.
- `protein`: DIA-NN protein group matrix, containing normalised quantities for 2476 protein groups in 434 single cells. Proteins are filtered at 1% FDR, using global q-values for protein groups and both global and run-specific q-values for precursors.

The `colData(brunner2022())` contains cell type annotations and batch annotations. The description of the `rowData` fields for the different assays can be found in the [DIA-NN documentation](#).

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation:** cells were detached with trypsin treatment, followed by strong pipetting, and isolate using FACS.
- **Sample preparation:** cell lysis by freeze-heat followed by sonication, overnight protein digestion with trypsin/lysC mix and desalting using EvoTips trap column (EvoSep)
- **Separation:** online EvoSep One LC system using a 5 cm x 75  $\mu\text{m}$  ID column with 1.9 $\mu\text{m}$  C18 beads (EvoSep) at 100nL/min flow rate.
- **Ionization:** 10 $\mu\text{m}$  ID zero dead volume electrospray emitter (Bruker Daltonik) + nanoelectrospray ion source (Captive spray, Bruker Daltonik)
- **Mass spectrometry:** DIA PASEF mode. Correlation between IM and m/z was used to synchronize the elution of precursors from each IM scan with the quadrupole isolation window. Five consecutive diaPASEF cycles. The collision energy was ramped linearly as a function of the IM from 59 eV at  $1/K0=1.6 \text{ Vs cm}^2$  to 20 eV at  $1/K0=0.6 \text{ Vs cm}^2$ .
- **Data analysis:** DIA-NN (1.8).

## Data collection

The data were collected from the PRIDE [repository](#) in the DIANN1.8\_SingleCells\_CellCycle.zip file.

We loaded the DIA-NN main report table and generated a sample annotation table based on the MS file names. We next combined the sample annotation and the DIANN tables into a `QFeatures::QFeatures` object following the scp data structure. We loaded the proteins group matrix as a `SingleCellExperiment::SingleCellExperiment` object, fixed ambiguous protein group names, and added the protein data as a new assay and link the precursors to proteins using the `Protein.Group` variable from the `rowData`.

## Source

The data were downloaded from PRIDE [repository](#) with accession ID PXD024043.

## References

Brunner, Andreas-David, Marvin Thielert, Catherine Vasilopoulou, Constantin Ammar, Fabian Coscia, Andreas Mund, Ole B. Hoerning, et al. 2022. "Ultra-High Sensitivity Mass Spectrometry Quantifies Single-Cell Proteome Changes upon Perturbation." *Molecular Systems Biology* 18 (3): e10798. [Link to article](#)

## Examples

brunner2022()

---

cong2020AC

*Cong et al. 2020 (Ana. Chem.): HeLa single cells*

---

## Description

Single-cell proteomics using the nanoPOTS sample processing device in combination with ultranarrow-bore (20um i.d.) packed-column LC separations and the Orbitrap Eclipse Tribrid MS. The dataset contains label-free quantitative information at PSM, peptide and protein level. The samples are single HeLa cells. Bulk samples (100 and 20 cells) were also included in the experiment to increase the identification rate thanks to between-run matching (cf MaxQuant).

## Usage

cong2020AC

## Format

A `QFeatures::QFeatures` object with 9 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- 100/20 HeLa cells: 2 assays containing PSM data for a bulk of 100 or 20 HeLa cells, respectively.
- Blank: assay containing the PSM data for a blank sample
- Single cell X: 4 assays containing PSM data for a single cell. The X indicates the replicate number.
- peptides: quantitative data for 12590 peptides in 7 samples (all runs combined).
- proteins: quantitative data for 1801 proteins in 7 samples (all runs combined).

Sample annotation is stored in `colData(cong2020AC())`.

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation:** The HeLa cells were diluted and aspirated using a microcapillary with a pulled tip.
- **Sample preparation** performed using the nanoPOTS device. Protein extraction using RapiGest (+ DTT) + alkylation (IAA) + Lys-C digestion + cleave RapiGest (formic acid)
- **Separation:** UltiMate 3000 RSLCnano pump with a home-packed nanoLC column (60cm x 20um i.d.; approx. 20 nL/min)
- **Ionization:** ESI (2,000V; Nanospray Flex)
- **Mass spectrometry:** Thermo Fisher Orbitrap Fusion Eclipse. MS1 settings: accumulation time = 246ms; resolution = 120,000; AGC = 1E6. MS/MS settings depend on quantity. All: AGC = 1E5. 20-100 cels: accumulation time = 246ms; resolution = 120,000. Single cells: accumulation time = 500ms; resolution = 240,000.
- **Data analysis:** MaxQuant (v1.6.3.3) + Excel

## Data collection

The PSM, peptide and protein data were collected from the PRIDE repository (accession ID: PXD016921). We downloaded the `evidence.txt` file containing the PSM identification and quantification results. The sample annotation was inferred from the samples names. The data were then converted to a `QFeatures::QFeatures` object using the `scp::readSCP()` function.

The peptide data were processed similarly from the `peptides.txt` file. The quantitative column names were adapted to match the PSM data. The peptide data were added to `QFeatures::QFeatures` object and link between the features were stored.

The protein data were similarly processed from the `proteinGroups.txt` file. The quantitative column names were adapted to match the PSM data. The peptide data were added to `QFeatures::QFeatures` object and link between the features were stored.

## Source

All files can be downloaded from the PRIDE repository PXD016921. The source link is: <ftp://ftp.pride.ebi.ac.uk/pride/data>

## References

Cong, Yongzheng, Yiran Liang, Khatereh Motamedchaboki, Romain Huguet, Thy Truong, Rui Zhao, Yufeng Shen, Daniel Lopez-Ferrer, Ying Zhu, and Ryan T. Kelly. 2020. "Improved Single-Cell Proteome Coverage Using Narrow-Bore Packed NanoLC Columns and Ultrasensitive Mass Spectrometry." *Analytical Chemistry*, January. ([link to article](#)).

## Examples

```
cong2020AC()
```

---

derks2022

*Derks et al. 2022 - plexDIA (Nat. Biotechnol.): PDAC vs melanoma cells vs monocytes*

---

## Description

Single cell proteomics data acquired by the Slavov Lab using the plexDIA protocol. It contains quantitative information from pancreatic ductal acinar cells (PDAC; HPAF-II), melanoma cells (WM989-A6-G3) and monocytes (U-937) at precursor and protein level. The each run acquired 3 samples thanks to mTRAQ multiplexing.

## Usage

```
derks2022
```

## Format

A `QFeatures::QFeatures` object with 66 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object. The assays either hold the DIA-NN main output report table or the DIA-NN MS1 extracted signal table. The DIA-NN main output report table contains the results of the spectrum identification and quantification. The DIA-NN MS1 extracted signal table contains quantification for

all mTRAQ channels if its precursors was identified in at least one of the channels, regardless of whether there is sufficient evidence in those channels at 1% FDR.

The data is composed of three datasets

1. **Bulk**: dataset containing bulk (100-cell) data acquired using a Q-Exactive mass spectrometer. Assays 1-3 contain data from the DIA-NN main output report; assay 4 is the DIA-NN MS1 extracted signal.
2. **tims**: dataset containing single-cell data acquired using a timsTOF-SCP mass spectrometer. Assays 5-15 contain data from the DIA-NN main output report; assay 16 is the DIA-NN MS1 extracted signal.
3. **qe**: dataset containing single-cell data acquired using a Q-Exactive mass spectrometer. Assays 17-64 contain data from the DIA-NN main output report; assay 65 is the DIA-NN MS1 extracted signal.

The last assay proteins contains the processed protein data table generated by the authors.

The `colData(derks2022())` contains cell type annotations and batch annotations. The description of the `rowData` fields for the different assays can be found in the [DIA-NN documentation](#).

### Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation**: CellenONE cell sorting.
- **Sample preparation** performed using the improved SCoPE2 protocol using the CellenONE liquid handling system. nPOP cell lysis (DMSO) + trypsin digestion + mTRAQ (3plex) labelling and pooling. A target library was generated as well to perform prioritized DDA (Huffman et al. 2022) using MaxQuant.Live (2.0.3).
- **Separation**: bulk - online nLC (Dionex UltiMate 3000 UHPLC) with a 25 cm × 75 μm IonOpticks Aurora Series UHPLC column (AUR2-25075C18A), 200nL/min. qe - online nLC (Dionex UltiMate 3000 UHPLC) with a 15 cm × 75 μm IonOpticks Aurora Series UHPLC column (AUR2-15075C18A), 200nL/min. tims - nanoElute liquid chromatography system (Bruker Daltonics) using a 25 cm × 75 μm, 1.6-μm C18 (AUR2-25075C18A-CSI, IonOpticks).
- **Ionization**: ESI.
- **Mass spectrometry**: cf article.
- **Data analysis**: DIA-NN (1.8.1 beta 16).

### Data collection

The data were collected from a shared Google Drive [folder](#) that is accessible from the SlavovLab website (see Source section).

For each dataset separately, we combined the sample annotation and the DIANN tables in a `QFeatures::QFeatures` object following the `scp` data structure. We then combined the three datasets in a single `QFeatures` object. We load the proteins table processed by the authors as a `SingleCellExperiment::SingleCellExperiment` object and adapted the sample names to match those in the `QFeatures` object. We added the protein data as a new assay and link the precursors to proteins using the `Protein.Group` variable from the `rowData`.

### Source

The data were downloaded from the [Slavov Lab](#) website. The raw data and the quantification data can also be found in the massIVE repository [MSV000089093](#).

## References

Derks, Jason, Andrew Leduc, Georg Wallmann, R. Gray Huffman, Matthew Willetts, Saad Khan, Harrison Specht, Markus Ralser, Vadim Demichev, and Nikolai Slavov. 2022. "Increasing the Throughput of Sensitive Proteomics by plexDIA." *Nature Biotechnology*, July. [Link to article](#)

## Examples

```
derks2022()
```

---

dou2019\_boosting

*Dou et al. 2019 (Anal. Chem.): testing boosting ratios*

---

## Description

Single-cell proteomics using nanoPOTS combined with TMT isobaric labeling. It contains quantitative information at PSM and protein level. The cell types are either "Raw" (macrophage cells), "C10" (epithelial cells), or "SVEC" (endothelial cells). Each cell is replicated 2 or 3 times. Each cell type was run using 3 levels of boosting: 0 ng (no boosting), 5 ng or 50 ng. When boosting was applied, 1 reference well and 1 boosting well were added, otherwise 1 empty well was added. Each boosting setting (0ng, 5ng, 50ng) was run in duplicate.

## Usage

```
dou2019_boosting
```

## Format

A `QFeatures::QFeatures` object with 7 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `Boosting_X_run_Y`: PSM data with 10 columns corresponding to the TMT-10plex channels. The X indicates the boosting amount (0ng, 5ng or 50ng) and Y indicates the run number (1 or 2).
- `peptides`: peptide data containing quantitative data for 13,462 peptides in 60 samples (run 1 and run 2 combined).
- `proteins`: protein data containing quantitative data for 1436 proteins and 60 samples (all runs combined).

Sample annotation is stored in `colData(dou2019_boosting())`.

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation:** single-cells from the three murine cell lines were isolated using FACS (BD Influx II cell sorter). Boosting sample were prepared (presumably in bulk) from 1:1:1 mix of the three cell lines.

- **Sample preparation** performed using the nanoPOTs device. Protein extraction (DMM + TCEAP) + alkylation (IAA) + Lys-C digestion + trypsin digestion + TMT-10plex labeling and pooling.
- **Separation:** nanoLC (Dionex UltiMate with an in-house packed 50cm x 30um LC columns; 50nL/min)
- **Ionization:** ESI (2,000V)
- **Mass spectrometry:** Thermo Fisher Orbitrap Fusion Lumos Tribrid (MS1 accumulation time = 50ms; MS1 resolution = 120,000; MS1 AGC = 1E6; MS2 accumulation time = 246ms; MS2 resolution = 60,000; MS2 AGC = 1E5)
- **Data analysis:** MS-GF+ + MASIC (v3.0.7111) + RomicsProcessor (custom R package)

### Data collection

The PSM data were collected from the MassIVE repository MSV000084110 (see Source section). The downloaded files are:

- Boosting\_\*ng\_run\_\*\_msgfplus.mzid: the MS-GF+ identification result files.
- Boosting\_\*ng\_run\_\*\_ReporterIons.txt: the MASIC quantification result files.

For each batch, the quantification and identification data were combined based on the scan number (common to both data sets). The combined datasets for the different runs were then concatenated feature-wise. To avoid data duplication due to ambiguous matching of spectra to peptides or ambiguous mapping of peptides to proteins, we combined ambiguous peptides to peptides groups and proteins to protein groups. Feature annotations that are not common within a peptide or protein group are separated by a ;. The sample annotation table was manually created based on the available information provided in the article. The data were then converted to a `QFeatures::QFeatures` object using the `scp::readSCP()` function.

We generated the peptide data. First, we removed PSM matched to contaminants or decoy peptides and ensured a 1% FDR. We aggregated the PSM to peptides based on the peptide (or peptide group) sequence(s) using the median PSM intensity. The peptide data for the different runs were then joined in a single assay (see `QFeatures::joinAssays`), again based on the peptide sequence(s). We then removed the peptide groups. Links between the peptide and the PSM data were created using `QFeatures::addAssayLink`. Note that links between PSM and peptide groups are not stored.

The protein data were downloaded from Supporting information section from the publisher's website (see Sources). The data is supplied as an Excel file `ac9b03349_si_004.xlsx`. The file contains 7 sheets from which we took the 2nd, 4th and 6th sheets (named 01 - No Boost raw data, 03 - 5ng boost raw data, 05 - 50ng boost raw data, respectively). The sheets contain the combined protein data for the duplicate runs given the boosting amount. We joined the data for all boosting ration based on the protein name and converted the data to a `SingleCellExperiment::SingleCellExperiment` object. We then added the object as a new assay in the `QFeatures::QFeatures` dataset (containing the PSM data). Links between the proteins and the corresponding PSM were created. Note that links to protein groups are not stored.

### Source

The PSM data can be downloaded from the massIVE repository MSV000084110. FTP link: <ftp://massive.ucsd.edu/MSV000084110>

The protein data can be downloaded from the [ACS Publications](#) website (Supporting information section).

## References

Dou, Maowei, Jeremy Clair, Chia-Feng Tsai, Kerui Xu, William B. Chrisler, Ryan L. Sontag, Rui Zhao, et al. 2019. “High-Throughput Single Cell Proteomics Enabled by Multiplex Isobaric Labeling in a Nanodroplet Sample Preparation Platform.” *Analytical Chemistry*, September ([link to article](#)).

## See Also

[dou2019\\_lysates](#), [dou2019\\_mouse](#)

## Examples

```
dou2019_boosting()
```

---

dou2019\_lysates

*Dou et al. 2019 (Anal. Chem.): HeLa lysates*

---

## Description

Single-cell proteomics using nanoPOTS combined with TMT multiplexing. It contains quantitative information at PSM and protein level. The samples are commercial HeLa lysates diluted to single-cell amounts (0.2 ng). The boosting wells contain the same digest but at higher amount (10 ng).

## Usage

```
dou2019_lysates
```

## Format

A `QFeatures::QFeatures` object with 3 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `HeLa_run_1`: PSM data with 10 columns corresponding to the TMT-10plex channels. Columns hold quantitative information for HeLa lysate samples (either 0, 0.2 or 10ng). This is the data for run 1.
- `HeLa_run_2`: PSM data with 10 columns corresponding to the TMT-10plex channels. Columns hold quantitative information for HeLa lysate samples (either 0, 0.2 or 10ng). This is the data for run 2.
- `peptides`: peptide data containing quantitative data for 13,934 peptides in 20 samples (run 1 and run 2 combined).
- `proteins`: protein data containing quantitative data for 1641 proteins in 20 samples (run 1 and run 2 combined).

Sample annotation is stored in `colData(dou2019_lysates())`.

### Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation:** commercially available HeLa protein digest (Thermo Scientific).
- **Sample preparation** performed using the nanoPOTs device. Protein extraction (DMM + TCEAP) + alkylation (IAA) + Lys-C digestion + trypsin digestion + TMT-10plex labeling and pooling.
- **Separation:** nanoLC (Dionex UltiMate with an in-house packed 50cm x 30um LC columns; 50nL/min)
- **Ionization:** ESI (2,000V)
- **Mass spectrometry:** Thermo Fisher Orbitrap Fusion Lumos Tribrid (MS1 accumulation time = 50ms; MS1 resolution = 120,000; MS1 AGC = 1E6; MS2 accumulation time = 246ms; MS2 resolution = 60,000; MS2 AGC = 1E5)
- **Data analysis:** MS-GF+ + MASIC (v3.0.7111) + RomicsProcessor (custom R package)

### Data collection

The PSM data were collected from the MassIVE repository MSV000084110 (see Source section). The downloaded files are:

- `HeLa_run*_msgfplus.mzid`: the MS-GF+ identification result files
- `HeLa_run*_ReporterIons.txt`: the MASIC quantification result files

For each batch, the quantification and identification data were combined based on the scan number (common to both data sets). The combined datasets for the different runs were then concatenated feature-wise. To avoid data duplication due to ambiguous matching of spectra to peptides or ambiguous mapping of peptides to proteins, we combined ambiguous peptides to peptides groups and proteins to protein groups. Feature annotations that are not common within a peptide or protein group are separated by a ;. The sample annotation table was manually created based on the available information provided in the article. The data were then converted to a `QFeatures::QFeatures` object using the `scp::readSCP()` function.

We generated the peptide data. First, we removed PSM matched to contaminants or decoy peptides and ensured a 1% FDR. We aggregated the PSM to peptides based on the peptide (or peptide group) sequence(s) using the median PSM intensity. The peptide data for the different runs were then joined in a single assay (see `QFeatures::joinAssays`), again based on the peptide sequence(s). We then removed the peptide groups. Links between the peptide and the PSM data were created using `QFeatures::addAssayLink`. Note that links between PSM and peptide groups are not stored.

The protein data were downloaded from Supporting information section from the publisher's website (see Sources). The data is supplied as an Excel file `ac9b03349_si_003.xlsx`. The file contains 7 sheets from which we only took the sheet 6 (named 5 - Run 1 and 2 raw data) with the combined protein data for the two runs. We converted the data to a `SingleCellExperiment::SingleCellExperiment` object and added the object as a new assay in the `QFeatures::QFeatures` dataset (containing the PSM data). Links between the proteins and the peptides were created. Note that links to protein groups are not stored.

### Source

The PSM data can be downloaded from the massIVE repository MSV000084110. FTP link: <ftp://massive.ucsd.edu/MSV>

The protein data can be downloaded from the [ACS Publications](#) website (Supporting information section).

## References

Dou, Maowei, Jeremy Clair, Chia-Feng Tsai, Kerui Xu, William B. Chrisler, Ryan L. Sontag, Rui Zhao, et al. 2019. "High-Throughput Single Cell Proteomics Enabled by Multiplex Isobaric Labeling in a Nanodroplet Sample Preparation Platform." *Analytical Chemistry*, September ([link to article](#)).

## See Also

[dou2019\\_mouse](#), [dou2019\\_boosting](#)

## Examples

```
dou2019_lysesates()
```

---

dou2019\_mouse

*Dou et al. 2019 (Anal. Chem.): murine cell lines*

---

## Description

Single-cell proteomics using nanoPOTS combined with TMT isobaric labeling. It contains quantitative information at PSM and protein level. The cell types are either "Raw" (macrophage cells), "C10" (epithelial cells), or "SVEC" (endothelial cells). Out of the 132 wells, 72 contain single cells, corresponding to 24 C10 cells, 24 RAW cells, and 24 SVEC. The other wells are either boosting channels (12), empty channels (36) or reference channels (12). Boosting and reference channels are balanced (1:1:1) mixes of C10, SVEC, and RAW samples at 5 ng and 0.2 ng, respectively. The different cell types were evenly distributed across 4 nanoPOTS chips. Samples were 11-plexed with TMT labeling.

## Usage

```
dou2019_mouse
```

## Format

A `QFeatures::QFeatures` object with 13 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `Single_Cell_Chip_X_Y`: PSM data with 11 columns corresponding to the TMT channels (see Notes). The X indicates the chip number (from 1 to 4) and Y indicates the row name on the chip (from A to C).
- `peptides`: peptide data containing quantitative data for 15,492 peptides in 132 samples (run 1 and run 2 combined).
- `proteins`: protein data containing quantitative data for 2331 proteins in 132 samples (all runs combined).

Sample annotation is stored in `colData(dou2019_mouse())`.

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation:** single-cells from the three murine cell lines were isolated using FACS (BD Influx II cell sorter).
- **Sample preparation** performed using the nanoPOTs device. Protein extraction (DMM + TCEAP) + alkylation (IAA) + Lys-C digestion + trypsin digestion + TMT-10plex labeling and pooling.
- **Separation:** nanoLC (Dionex UltiMate with an in-house packed 50cm x 30um LC columns; 50nL/min)
- **Ionization:** ESI (2,000V)
- **Mass spectrometry:** Thermo Fisher Orbitrap Fusion Lumos Tribrid (MS1 accumulation time = 50ms; MS1 resolution = 120,000; MS1 AGC = 1E6; MS2 accumulation time = 246ms; MS2 resolution = 60,000; MS2 AGC = 1E5)
- **Data analysis:** MS-GF+ + MASIC (v3.0.7111) + RomicsProcessor (custom R package)

## Data collection

The PSM data were collected from the MassIVE repository MSV000084110 (see Source section). The downloaded files are:

- `Single_Cell_Chip_*_*_msgfplus.mzid`: the MS-GF+ identification result files.
- `Single_Cell_Chip_*_*_ReporterIons.txt`: the MASIC quantification result files.

For each batch, the quantification and identification data were combined based on the scan number (common to both data sets). The combined datasets for the different runs were then concatenated feature-wise. To avoid data duplication due to ambiguous matching of spectra to peptides or ambiguous mapping of peptides to proteins, we combined ambiguous peptides to peptides groups and proteins to protein groups. Feature annotations that are not common within a peptide or protein group are separated by a ;. The sample annotation table was manually created based on the available information provided in the article. The data were then converted to a `QFeatures::QFeatures` object using the `scp::readSCP()` function.

We generated the peptide data. First, we removed PSM matched to contaminants or decoy peptides and ensured a 1% FDR. We aggregated the PSM to peptides based on the peptide (or peptide group) sequence(s) using the median PSM intensity. The peptide data for the different runs were then joined in a single assay (see `QFeatures::joinAssays`), again based on the peptide sequence(s). We then removed the peptide groups. Links between the peptide and the PSM data were created using `QFeatures::addAssayLink`. Note that links between PSM and peptide groups are not stored.

The protein data were downloaded from Supporting information section from the publisher's website (see Sources). The data is supplied as an Excel file `ac9b03349_si_005.xlsx`. The file contains 7 sheets from which we only took the 2nd (named 01 - Raw sc protein data) with the combined protein data for the 12 runs. We converted the data to a `SingleCellExperiment::SingleCellExperiment` object and added the object as a new assay in the `QFeatures::QFeatures` dataset (containing the PSM data). Links between the proteins and the corresponding PSM were created. Note that links to protein groups are not stored.

## Note

Although a TMT-10plex labeling is reported in the article, the PSM data contained 11 channels for each run. Those 11th channel contain mostly missing data and are hence assumed to be empty channels.

## Source

The PSM data can be downloaded from the massIVE repository MSV000084110. FTP link: <ftp://massive.ucsd.edu/MSV>

The protein data can be downloaded from the [ACS Publications](#) website (Supporting information section).

## References

Dou, Maowei, Jeremy Clair, Chia-Feng Tsai, Kerui Xu, William B. Chrisler, Ryan L. Sontag, Rui Zhao, et al. 2019. “High-Throughput Single Cell Proteomics Enabled by Multiplex Isobaric Labeling in a Nanodroplet Sample Preparation Platform.” *Analytical Chemistry*, September ([link to article](#)).

## See Also

[dou2019\\_lysates](#), [dou2019\\_boosting](#)

## Examples

```
dou2019_mouse()
```

---

gregoire2023\_mixCTRL *Grégoire et al. 2023 - mixCTRL (arXiv): benchmark using monocytes/macrophages*

---

## Description

Single cell proteomics data acquired using the SCoPE2 protocol. The dataset contains two monocytes cell lines (THP1 and U937) as well as controlled mixtures of both and macrophage-like cells produced upon PMA treatment. It contains quantitative information at PSM, peptide and protein levels. Data was acquired using Lumos Orbitrap (mainly) and timsTOF SCP mass spectrometers.

## Usage

```
gregoire2023_mixCTRL
```

## Format

A `QFeatures::QFeatures` object with 119 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- Assays 1-42: PSM data acquired with a TMT-16plex protocol, hence those assays contain 16 columns. Columns hold quantitative information from single-cell channels, carrier channels, blank (negative control) channels and unused channels.
- Assays 43-84: peptide data resulting from the PSM to peptide aggregation of the 42 PSM assays.
- Assays 85-91: peptide data for each of the 7 acquisition batches. Peptide data were joined based on their respective acquisition batches.
- Assays 92-98: normalised peptide data.

- Assays 99-105: normalised and log-transformed peptide data.
- Assays 106-112: protein data for each of the 7 acquisition batches. Normalised and log-transformed peptide data were aggregated to protein.
- Assays 113-119: Batch corrected protein data. Normalised and log-transformed protein data were batch corrected to remove technical variability induced by runs and channels.

All the data has been filtered to keep high quality features and samples.

The `colData(gregoire2023_mixCTRL())` contains cell type annotation and batch annotation that are common to all assays. The description of the `rowData` fields for the PSM data can be found in the [sage documentation](#).

### Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see *References*).

- **Cell isolation:** BD FACSAria III cell sorting.
- **Sample preparation** performed using the SCoPE2 protocol: mPOP cell lysis + trypsin digestion + TMT-16plex labeling and pooling.
- **Separation:** online nLC (Ultimate 3000 LC System or Vanquish Neo UHPLC System) with a BioZen Peptide Polar C18 250 x 0.0075mm column.
- **Mass spectrometry:** Orbitrap Fusion Lumos Tribrid (MS1 resolution = 70,000; MS2 accumulation time = 120ms; MS2 resolution = 70,000) and timsTOF SCP.
- **Data preprocessing:** Sage.

### Data collection

The PSM data were collected from a Zenodo archive (see Source section). The folder contains the following files of interest:

- `results.sage.cbio.tsv`: the sage identification output file for batches acquired on the Lumos MS.
- `results.sage.giga.tsv`: the sage identification output file for batches acquired on the timsTOF SCP MS.
- `quant.cbio.tsv`: the sage quantification output file for batches acquired on the Lumos MS.
- `quant.giga.tsv`: the sage quantification output file for batches acquired on the timsTOF SCP MS.
- `sampleAnnotation_batch.csv`: sample annotation for each acquisition batch. There are in total 8 different annotation files.

We combined the sample annotations in a single table. We also combined `cbio` and `giga` tables together and merged resulting identification and quantification tables. Both annotation and features tables are then combined in a single `QFeatures::QFeatures` object using the `scp::readSCP()` function.

The `QFeatures::QFeatures` object was processed as described in the author's manuscript (see source). Note that the imputed assays were used in the paper for illustrative purposes only and have not been reproduced here.

### Source

The data were downloaded from the [Zenodo repository](#). The raw data and the quantification data can also be found in the ProteomeXchange Consortium via the [PRIDE partner repository](#), project PXD046211.

## References

Samuel Grégoire, Christophe Vanderaa, Sébastien Pyr dit Ruys, Gabriel Mazzucchelli, Christopher Kune, Didier Vertommen and Laurent Gatto. 2023. *Standardised workflow for mass spectrometry-based single-cell proteomics data processing and analysis using the scp package*. arXiv. DOI:[10.48550/arXiv.2310.13598](https://doi.org/10.48550/arXiv.2310.13598)

## Examples

```
gregoire2023_mixCTRL()
```

---

guise2024

*Guise et al. 2020 (Cell Rep.): postmortem ALS spinal moto neurons*

---

## Description

Single-cell proteomics data from postmortem human spinal moto neurons (MN) obtained from control donors or donors with amyotrophic lateral sclerosis (ALS). The data were generated following the NanoPOTS protocol. Cells were isolated from samples obtained by the university of Miami Brain Bank using laser capture microdissection (LCM). Additional information about the amount of TDP-43 intra-cellular levels has been assigned into levels 0 to 4.

## Usage

```
guise2024
```

## Format

A `QFeatures::QFeatures` object with 102 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `F*`: 100 assays containing PSM data.
- `peptides`: quantitative data for 34,315 peptides in 108 samples. All samples combined, along with 8 additional unannotated samples.
- `proteins`: quantitative data for 4,437 protein groups in 108 samples. All samples combined, along with 8 additional unannotated samples.

Sample annotation is stored in `colData(guise2024())`.

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation**: The MN were isolated from samples obtained by the university of Miami Brain Bank using LCM.
- **Sample preparation** performed using the nanoPOTS workflow. Cells are treated with 0.1% DDM (for lysis) added with DTT (protein reduction), then IAA (alkylation), then Lys-C and trypsin (protein digestion).

- **Separation:** Samples were injected on the column using an Ultimate 3000 RSLCnano pump. The in-line loading column is a home-packed SPE column (5cm x 75um) while the peptide separation is performed on a an in-house-packed analytical SPE column (50 cm x 30um), using a 20nL/min flow rate.
- **Ionization:** nanospray emitter (2,000V)
- **Mass spectrometry:** Orbitrap Exploris 480. HCD fragmentation. MS1 settings: accumulation time = 200 ms; resolution = 120,000; AGC = 1E6. MS2 settings: exclusion duration = 90 s; accumulation time = 500 ms; resolution = 30,000; AGC = 1E5.
- **Data analysis:** Sequest HT in Proteome Discoverer (v2.5) and the search database is Swiss-Prot (July 2020).

## Data collection

All data were collected from the MassIVE repository (accession ID: MSV000092119).

The sample annotations were combined from the tables in Biogen\_TDP43\_Round2\_Reanalysis\_10-13-2021\_InputFiles and in Groups.txt.

The PSM data were found in the Biogen\_TDP43\_Round2\_Reanalysis\_10-13-2021\_PSMs.txt file. The data were converted to a `QFeatures::QFeatures` object using the `scp::readSCP()` function. We could not find sample annotations for MS run ID: F61, F34, F42, F88, F77, F8, F21, F5.

The peptide data were found in the Biogen\_TDP43\_Round2\_Reanalysis\_10-13-2021\_PeptideGroups.txt file. The column names holding the quantitative data were adapted to match the sample names in the `QFeatures::QFeatures` object. The data were then converted to a `SingleCellExperiment::SingleCellExperiment` object and then inserted in the `QFeatures::QFeatures` object.

A similar procedure was applied to the protein data. The data were found in the Biogen\_TDP43\_Round2\_Reanalysis\_10-13-2021\_ProteinGroups.txt file. The column names were adapted, the data were converted to a `SingleCellExperiment::SingleCellExperiment` object and then inserted in the `QFeatures::QFeatures` object.

## Source

All data can be downloaded from the MassIVE repository MSV000092119. The source link is: <ftp://massive.ucsd.edu/v05/MSV000092119/>

## References

Guise, Amanda J., Santosh A. Misal, Richard Carson, Jen-Hwa Chu, Hannah Boekweg, Daisha Van Der Watt, Nora C. Welsh, et al. 2024. "TDP-43-Stratified Single-Cell Proteomics of Postmortem Human Spinal Motor Neurons Reveals Protein Dynamics in Amyotrophic Lateral Sclerosis." Cell Reports 43 (1): 113636. ([link to article](#)).

## Examples

```
guise2024()
```

---

hu2023\_K562

*Hu et al, 2023 (The Journal of Physical Chemistry B): Correlated protein modules*

---

### Description

They demonstrate the correlations between the levels of pairs of proteins in single-cell proteomics (SCP) at steady state. In measuring pairwise correlations among 1000 proteins in a population of K562 cells and oocytes, they observed many correlated protein modules (CPMs) that are functionally involved in certain biological functions. Certain CPMs are specific to a particular cell type, some common to different cell types. Additionally, compared to single-cell transcriptomics and bulk proteomics, protein correlations are functionally and experimentally more significant in SCP than those corresponding mRNAs.

### Usage

hu2023\_K562

### Format

Two `SingleCellExperiment::SingleCellExperiment` objects:

- `proteins_K562`: protein data containing quantitative data for 1249 proteins and 69 single-cells with zero imputation.
- `proteins_oocyte`: protein data containing quantitative data for 3422 proteins and 137 single-cells with zero imputation.

The `colData(hu2023_oocyte())` contains cell type annotation. The `colData(hu2023_K562())` contains cell type annotation.

### Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation:** K562 cells were re-suspended and washed in cold PBS. Single cells/10 cells were sorted into 96-well plates using a FACSARIA instrument. Oocyte-cumulus complexes from C57/6J mice were collected after PMSG and HCG injections, with hyaluronidase used to remove cumulus cells. All samples stored at -80 degrees Celsius.
- **Sample preparation** Cells were digested with trypsin at 37 degrees Celsius for 3 hours. For label-free proteomics, digestion was terminated by adding 0.43% TFA and 1% ACN in water, followed by drying in a concentrator. Peptides were resuspended in 0.1% TFA and 1% ACN, and then transferred to sample tubes for LC-MS/MS analysis.
- **Separation:** 4 microliters of peptide digests were injected into a high-performance chromatography column (IonOpticks) and separated at a flow rate of 100 nL/min using a nanoflow liquid chromatography system. The effective gradient was 70 mins, allowing 16 cells per day.
- **Ionization:** Peptides were analyzed using an Orbitrap Eclipse mass spectrometer with a FAIMS Pro interface. FAIMS compensation voltages of -55 and -70 V were applied, with a 1-second cycle time for both voltages.

- **Mass spectrometry:** MS spectra were acquired with the Orbitrap analyzer, while MS/MS spectra were acquired with a linear ion trap analyzer. The maximum ion injection time for MS/MS was 200 ms.
- **Data analysis:** MS raw files were searched against the UniProt human protein database and an in-house contamination database using Proteome Discoverer(2.4). Label-free quantification was based on peak intensity with the match-between-runs (MBR) feature enabled.

## Data collection

The oocyte protein data shared by the author and it is accessible from the [Shared File](#) The K563 protein data is accessible from the GitHub (<https://github.com/dionezhang/CPM/blob/master/ProteinAbundance.Rdata>).

- DataMatrix-oocyte-20240614.csv: normalized imputed protein matrix
- ProteinAbundance.Rdata: protein matrices (normalized, log transformed)

We initialized an empty QFeatures object and added the corresponding protein assays as `SingleCellExperiment::SingleCellExperiment` objects.

The oocyte protein data were exported from the shared link as (DataMatrix-oocyte-20240614.csv). The data were formatted to a `SingleCellExperiment::SingleCellExperiment` object and the SampleType information were added as only metadata, and stored in the colData. The object is then added to the `QFeatures::QFeatures` object.

The 562 cells protein data were downloaded from the GitHub link and loaded to the memory. The Norm object were formatted to a `SingleCellExperiment::SingleCellExperiment` object and the SampleType information were added as only metadata, and stored in the colData. The object is then added to the `QFeatures::QFeatures` object.

## Source

The oocyte data were downloaded from the [Shared File](#) The K563 cells protein data downloaded from the GitHub (<https://github.com/dionezhang/CPM/blob/master/ProteinAbundance.Rdata>) The raw data and the quantification data can also be found in the MassIVE repository MSV000089625: <ftp://MSV000089625@massive.ucsd.edu/>.

## References

- Hu, M., Zhang, Y., Yuan, Y., Ma, W., Zheng, Y., Gu, Q., & Xie, X. S. 2023. “Correlated protein modules revealing functional coordination of interacting proteins are detected by single-cell proteomics.”. The Journal of Physical Chemistry B, ([link to article](#)).

## Examples

```
hu2023_oocyte()  
hu2023_K562()
```

---

khan2023

*Khan et al, 2023 (biorRxiv): Epithelial–Mesenchymal Transition*

---

## Description

Single-cell samples were prepared using the nPOP sample preparation method. Proteomics data were acquired using the SCoPE2 protocol on a Thermo Scientific Q-Exactive mass spectrometer. The dataset contains quantitative information on 421 MCF-10A single cells undergoing epithelial–mesenchymal transition (EMT) triggered by TGF beta. The data are available at the PSM, and protein levels. The paper investigates the dynamics of correlation modules at the protein level.

## Usage

khan2023

## Format

A `QFeatures::QFeatures` object with 47 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- Assay 1-44: PSM data acquired with a TMTPro 16plex protocol, hence those assays contain 16 columns. Columns hold quantitative information from single-cell channels, carrier channels, reference channels, empty (negative control) channels and unused channels.
- peptides: peptide data containing quantitative data for 10055 peptides and 421 single-cells.
- proteins\_imputed: protein data containing quantitative data for 4571 proteins and 420 single-cells with k-nearest neighbors (KNN) imputation.
- proteins\_unimputed: protein data containing quantitative data for 4571 proteins and 420 single-cells without imputation.

The `colData(khan2023())` contains cell type and batch annotations that are common to all assays. The description of the `rowData` fields for the PSM data can be found in the [MaxQuant documentation](#).

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation:** CellenONE cell sorting.
- **Sample preparation** performed using the SCoPE2 protocol. nPOP cell lysis (DMSO) + trypsin digestion + TMTPro 16plex protocol.
- **Separation:** online nLC (DionexUltiMate 3000 UHPLC with a 25cm x 75um IonOpticks Odyssey Series column (ODY3-25075C18); 200nL/min).
- **Ionization:** ESI (1,700 V).
- **Mass spectrometry:** Thermo Scientific Q-Exactive (MS1 resolution = 70,000; MS1 accumulation time = 300ms; MS2 resolution = 70,000).
- **Data analysis:** MaxQuant(2.4.13.0) + DART-ID.

## Data collection

The PSM data were collected from a shared Google Drive folder that is accessible from the SlavovLab website (see Source section). The folder ('/002-singleCellDataGeneration') contains the following files of interest:

- `ev_updated_NS.DIA.txt`: the MaxQuant/DART-ID output file
- `annotation.csv`: sample annotation
- `batch.csv`: batch annotation

We combined the sample annotation and the batch annotation in a single table. We also formatted the quantification table so that columns match with those of the annotation and filter only for single-cell runs. Both table are then combined in a single `QFeatures::QFeatures` object using the `scp::readSCP()` function.

The peptide data were taken from the same google drive folder (`EpiToMesen.TGFB.nPoP_trial11_pepByCellMatrix_NS`). The data were formatted to a `SingleCellExperiment::SingleCellExperiment` object and the sample metadata were matched to the column names (mapping is retrieved after running the `SCoPE2 R` script, `EMTTGFB_singleCellProcessing.R`) and stored in the `colData`. The object is then added to the `QFeatures::QFeatures` object and the rows of the PSM data are linked to the rows of the peptide data based on the peptide sequence information through an `AssayLink` object.

The imputed protein data were taken from the same google drive folder (`EpiToMesen.TGFB.nPoP_trial11_1PercDartFD`). The data were formatted to a `SingleCellExperiment::SingleCellExperiment` object and the sample metadata were matched to the column names (mapping is retrieved after running the `SCoPE2 R` script, `EMTTGFB_singleCellProcessing.R`) and stored in the `colData`. The object is then added to the `QFeatures::QFeatures` object and the rows of the peptide data are linked to the rows of the protein data based on the protein sequence information through an `AssayLink` object.

The unimputed protein data were taken from the same google drive folder (`EpiToMesen.TGFB.nPoP_trial11_1PercDart`). The data were formatted and added exactly as imputed data.

## Source

The data were downloaded from the [Slavov Lab](#) website via a shared Google Drive [folder](#). The raw data and the quantification data can also be found in the MassIVE repository `MSV000092872`: `ftp://MSV000092872@massive.ucsd.edu/`.

## References

Saad Khan, Rachel Conover, Anand R. Asthagiri, Nikolai Slavov. 2023. "Dynamics of single-cell protein covariation during epithelial–mesenchymal transition." *Journal of Proteome Research*. ([link to article](#)).

## Examples

`khan2023()`

---

krull12024

*Krull et al, 2024 (Nature Communications): IFN- $\gamma$  response*

---

## Description

They develop a new strategy for data-independent acquisition (DIA) that leverages the co-analysis of low-input samples alongside a corresponding enhancer (ME) of higher input. Using DIA-ME, they investigate the proteomic response of U-2 OS cells to interferon gamma (IFN- $\gamma$ ) at the single-cell level.

## Usage

krull12024

## Format

A `QFeatures::QFeatures` object with 159 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object.

- Assay 1-158: DIA-NN main output report table split for each acquisition run. First 15 run acquires 10 single cells (MEs) and, remaining 143 run acquires 1 single cell. It contains the results of the spectrum identification and quantification.
- proteins: DIA-NN protein group matrix, containing normalised quantities for 1553 protein groups in 143 single cells. Proteins are filtered at (`Q.Value <= 0.01`), (`Lib.Q.Value <= 0.01`), and (`Lib.PG.Q.Value <= 0.01`).

The `colData(krull12024())` contains cell type annotations. The description of the `rowData` fields for the different assays can be found in the [DIA-NN documentation](#).

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation:** cells were detached with trypsin digestion, followed by dilution in 1.5 mL PBS, and isolated using BD FACSAria III instrument.
- **Sample preparation:** Sorted single cells were collected in lysis buffer (50 mM TEAB, pH 8.5, and 0.025% DDM), denatured at 70 degrees Celsius for 30 minutes. Samples were acidified with 0.5% FA and transferred to auto sampler plates for mass spectrometry analysis.
- **Separation:** Peptides were injected in a 2 microliter volume onto a (25 cm x 75 micrometer) ID column at a flow rate of 300 nL/min, separated using a gradient of ACN in water with 0.1% FA over 15 minutes, connected to a nano-ESI source.
- **Ionization:** Ionization was performed using a 1,500 V capillary voltage with 3.0 L/min dry gas and a dry temperature of 180 degrees Celsius. MS data acquisition was conducted in diaPASEF mode using a timsTOF Pro mass spectrometer.
- **Mass spectrometry:** MS1 scans covered a range of 200-1,700 m/z, while DIA window isolation targeted 475-1,000 m/z with eight DIA scans per cycle. Fragmentation was triggered by collision energy ranging from 45 eV to 27 eV depending on the ion mobility.
- **Data analysis:** Data was processed using DIA-NN (v1.8.0) and Spectronaut 18 in a library-free approach, using deep learning for spectrum prediction, retention times, and ion mobility.

## Data collection

The data were collected from the PRIDE [repository](#) in the 03\_SingleCell\_Searches.zip file.

We loaded the DIA-NN main report table and generated a sample annotation table based on the MS file names. We next combined the sample annotation and the DIANN tables into a `QFeatures::QFeatures` object following the scp data structure. We loaded the proteins group matrix as a `SingleCellExperiment::SingleCellExperiment` object, and added the protein data as a new assay and link the precursors to proteins using the `Protein.Group` variable from the `rowData`.

## Source

The data were downloaded from PRIDE [repository](#) with accession ID PXD053464.

## References

Krull, K. K., Ali, S. A., & Krijgsveld, J. 2024. "Enhanced feature matching in single-cell proteomics characterizes IFN- $\gamma$  response and co-existence of Cell States." *Nature Communications*, 15(1). [Link to article](#)

## Examples

```
krull12024()
```

---

leduc2022

*Deprecated leduc2022 dataset*

---

## Description

The leduc2022 dataset has been updated to include plexDIA and pSCoPE data. The new datasets names are leduc2022\_pSCoPE (previously leduc2022) and leduc2022\_plexDIA (new). See the respective documentation pages for more information.

## Usage

```
leduc2022()
```

## Value

The leduc2022\_pSCoPE dataset.

---

`leduc2022_plexDIA`*Leduc et al. 2022 - plexDIA (biorRxiv): melanoma cells*

---

## Description

Single cell proteomics data acquired by the Slavov Lab. This is the dataset associated to the fourth version of the preprint (and the Genome Biology publication). It contains quantitative information of melanoma cells at precursor, peptide and protein level. This version of the data was acquired using the plexDIA MS acquisition protocol.

## Usage

`leduc2022_plexDIA`

## Format

A `QFeatures::QFeatures` object with 48 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- **Assay 1-45:** precursor data acquired with a mTRAQ-3 protocol, hence those assays contain 3 columns. Columns hold quantitative information from single cells or negative control samples.
- **Ms1Extracted:** the DIA-NN MS1 extracted signal, it combines the information from assays 1-45.
- **peptides:** peptide data containing quantitative data for 3,608 peptides and 104 single cells. The data were filtered to 1% protein FDR.
- **proteins:** protein data containing quantitative data for 508 proteins and 105 single cells. Note that the peptide and protein data provided by the authors differ by 3 samples. The precursor data were aggregated to protein intensity using maxLFQ. The protein data were further median normalized by column and by row, log2 transformed, impute using KNN (k = 3), again median normalized by column and by row, batch corrected using ComBat, and median normalized by column and by row once more.

The `colData(leduc2022_plexDIA())` contains cell type annotation and batch annotation that are common to all assays. The description of the `rowData` fields for the precursor data can be found in the [DIA-NN documentation](#).

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation:** CellenONE cell sorting.
- **Sample preparation** performed using the improved SCoPE2 protocol using the CellenONE liquid handling system. nPOP cell lysis (DMSO) + trypsin digestion + mTRAQ-3 labeling and pooling.
- **Separation:** online nLC (DionexUltiMate 3000 UHPLC with a 25cm x 75um IonOpticks Aurora Series UHPLC column; 200nL/min).
- **Ionization:** ESI (1,800V).

- **Mass spectrometry:** Thermo Scientific Q-Exactive. The duty cycle = 1 MS1 + 4 DIA MS2 windows (120 Th, 120 Th, 200 Th and 580 Th, spanning 378-1,402 m/z). Each MS1 and MS2 scan was conducted at 70,000 resolving power,  $3 \times 10^6$  AGC and 300ms maximum injection time.
- **Data analysis:** DIA-NN.

### Data collection

The PSM data were collected from a shared Google Drive folder that is accessible from the SlavovLab website (see Source section). The folder contains the following files of interest:

- `annotation_plexDIA.csv`: sample annotation
- `report_plexDIA_me1_nPOP.tsv`: the DIA-NN output file with the precursor data
- `report.pr_matrix_channels_ms1_extracted.tsv`: the DIA-NN output file with the combined precursor data
- `plexDIA_peptide.csv`: the processed data table containing the peptide data
- `plexDIA_protein_imputed.csv`: the processed data table containing the protein data

We removed the failed runs as identified by the authors. We also formatted the annotation and precursor quantification tables to facilitate matching between corresponding columns. Both annotation and quantification tables are then combined in a single `QFeatures::QFeatures` object using `scp::readSCPfromDIANN()`.

The `plexDIA_peptide.csv` and `plexDIA_protein_imputed.csv` files were loaded and formatted as `SingleCellExperiment::SingleCellExperiment` objects. The column names were adapted to match those in the `QFeatures` object. The `SingleCellExperiment` objects were then added to the `QFeatures::QFeatures` object and the rows of the peptide data are linked to the rows of the precursor data based on the peptide sequence or the protein name through an `AssayLink` object.

### Source

The links to the data were found on the [Slavov Lab website](#). The data were downloaded from the [Google drive folder 1](#) and [Google drive folder 2](#). The raw data and the quantification data can also be found in the massIVE repository MSV000089159: <ftp://massive.ucsd.edu/MSV000089159>.

### References

Andrew Leduc, Gray Huffman, and Nikolai Slavov. 2022. “Droplet Sample Preparation for Single-Cell Proteomics Applied to the Cell Cycle.” *bioRxiv*. [Link to article](#)

Andrew Leduc, Gray Huffman, Joshua Cantlon, Saad Khan, and Nikolai Slavov. 2022. “Exploring Functional Protein Covariation across Single Cells Using nPOP.” *Genome Biology* 23 (1): 261. [Link to article](#)

Jason Derks, Andrew Leduc, Georg Wallmann, Gray Huffman, Matthew Willetts, Saad Khan, Harrison Specht, Markus Ralser, Vadim Demichev, and Nikolai Slavov. 2023. “Increasing the Throughput of Sensitive Proteomics by plexDIA.” *Nature Biotechnology* 41 (1): 50–59. [Link to article](#)

### See Also

[leduc2022\\_pSCoPE](#)

**Examples**

```
leduc2022_plexDIA()
```

---

```
leduc2022_pSCoPE
```

```
Leduc et al. 2022 - pSCoPE (biorRxiv): melanoma cells vs monocytes
```

---

**Description**

Single cell proteomics data acquired by the Slavov Lab. This is the dataset associated to the third version of the preprint. It contains quantitative information of melanoma cells and monocytes at PSM, peptide and protein level. This version of the data was acquired using the pSCoPE MS acquisition approach.

**Usage**

```
leduc2022_pSCoPE
```

**Format**

A `QFeatures::QFeatures` object with 138 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- Assay 1-134: PSM data acquired with a TMT-18plex protocol, hence those assays contain 18 columns. Columns hold quantitative information from single-cell channels, carrier channels, reference channels, empty (negative control) channels and unused channels.
- `peptides`: peptide data containing quantitative data for 20,804 peptides and 1556 single-cells. These data have been filtered to keep high-quality PSMs, all batches have been normalized to the reference channel, PSMs were aggregated to peptides, and single-cells with low median coefficient of variation were kept.
- `peptides_log`: peptide data containing quantitative data for 12,284 peptides and 1543 single-cells. The peptides data was further normalized, highly missing peptides were removed and the quantifications were log-transformed.
- `proteins_norm2`: protein data containing quantitative data for 2844 proteins and 1543 single-cells. The peptides from `peptides_log` were aggregated to proteins and normalized.
- `proteins_processed`: protein data containing quantitative data for 2844 proteins and 1543 single-cells. The `proteins_norm2` data were imputed, batch corrected and normalized.

The `colData(leduc2022_pSCoPE())` contains cell type annotation, LC batch information, the TMT label, the MS run ID. We also added the sample prep annotations provided by the cellenONE dispensing device (only for single cells): time stamp of cell isolation by the device, the diameter and elongation of the cell, the ID of the sample glass side (4 slides in total), the field within the glass (each slide is divided in 4 field), the pooled well ID (each field contains 9 pools), the x and y coordinates of each cell dropped in a field and of each cell pool upon pickup. Finally, we also retrieved the melanoma subpopulation generated by the authors upon data analysis. The main population is encoded as A while the small population is encoded B. The description of the `rowData` fields for the PSM data can be found in the [MaxQuant documentation](#).

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation:** CellenONE cell sorting.
- **Sample preparation** performed using the improved SCoPE2 protocol using the CellenONE liquid handling system. nPOP cell lysis (DMSO) + trypsin digestion + TMT-18plex labeling and pooling. A target library was generated as well to perform prioritized DDA (Huffman et al. 2022) using MaxQuant.Live (2.0.3).
- **Separation:** online nLC (DionexUltiMate 3000 UHPLC with a 25cm x 75um IonOpticks Aurora Series UHPLC column; 200nL/min).
- **Ionization:** ESI (1,800V).
- **Mass spectrometry:** Thermo Scientific Q-Exactive (MS1 resolution = 70,000; MS2 accumulation time = 300ms; MS2 resolution = 70,000). Prioritized data acquisition was performed using the pSCoPE protocol (Huffman et al. 2022)
- **Data analysis:** MaxQuant (1.6.17.0) + DART-ID

## Data collection

The PSM data were collected from a shared Google Drive folder that is accessible from the SlavovLab website (see Source section). The folder contains the following files of interest:

- `ev_updated.txt`: the MaxQuant/DART-ID output file
- `annotation.csv`: sample annotation
- `batch.csv`: batch annotation
- `t0.csv`: the processed data table containing the peptides data
- `t3.csv`: the processed data table containing the peptides\_log data
- `t4b.csv`: the processed data table containing the proteins\_norm2 data
- `t6.csv`: the processed data table containing the proteins\_processed data

We combined the sample annotation and the batch annotation in a single table. We also formatted the quantification table so that columns match with those of the annotations. Both annotation and quantification tables are then combined in a single `QFeatures::QFeatures` object using the `scp::readSCP()` function.

The 4 CSV files were loaded and formatted as `SingleCellExperiment::SingleCellExperiment` objects and the sample metadata were matched to the column names (mapping is retrieved after running the author's original R script) and stored in the `colData`. The object is then added to the `QFeatures::QFeatures` object (containing the PSM assays) and the rows of the peptide data are linked to the rows of the PSM data based on the peptide sequence information through an `AssayLink` object.

## Source

The data were downloaded from the [Slavov Lab](#) website. The raw data and the quantification data can also be found in the massIVE repository MSV000089159: <ftp://massive.ucsd.edu/MSV000089159>.

## References

- Andrew Leduc, Gray Huffman, and Nikolai Slavov. 2022. "Droplet Sample Preparation for Single-Cell Proteomics Applied to the Cell Cycle." *bioRxiv*. [Link to article](#)
- Gray Huffman, Andrew Leduc, Christoph Wichmann, Marco di Gioia, Francesco Borriello, Harrison Specht, Jason Derks, et al. 2022. "Prioritized Single-Cell Proteomics Reveals Molecular and Functional Polarization across Primary Macrophages." *bioRxiv*. [Link to article](#).

**See Also**[leduc2022\\_plexDIA](#)**Examples**`leduc2022_pSCoPE()`

---

`liang2020_hela`*Liang et al. 2020 (Anal. Chem.): HeLa cells (MaxQuant preprocessing)*

---

**Description**

Single-cell proteomics data from HeLa cells using the autoPOTS acquisition workflow. The samples contain either no cells (blanks), 1 cell, 10 cells, 150 cells or 500 cells. Samples containing between 0 and 10 cells are isolated using micro-pipetting while samples containing between 150 and 500 cells were prepared using dilution of a bulk sample.

**Usage**`liang2020_hela`**Format**

A `QFeatures::QFeatures` object with 17 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `HeLa_*`: 15 assays containing PSM data.
- `peptides`: quantitative data for 48705 peptides in 15 samples (all runs are combined).
- `proteins`: quantitative data for 3970 protein groups in 15 samples (all runs combined).

Sample annotation is stored in `colData(liang2020_hela())`.

**Acquisition protocol**

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation:** The HeLa cells come from a commercially available cell line. Samples containing between 0 and 10 cells were isolated using micro-manipulation and the counts were validated using a microscope. Samples containing between 150 and 500 cells were prepared by diluting a bulk sample and the exact counts were evaluated by obtaining photomicrographs.
- **Sample preparation** performed using the autoPOTS workflow that relied on the OT-2 pipetting robot. Cells are lysed using sonication. Samples are then processed by successive incubation with DTT (reduction), then IAA (alkylation), then Lys-C and trypsin (protein digestion).
- **Separation:** Samples were injected on the column using a modified Ultimate WPS-3000 TPL autosampler coupled to an UltiMate 3000 RSLCnano pump. The LC column is a home-packed nanoLC column (45cm x 30µm; 40nL/min)
- **Ionization:** Nanospray Flex ion source (2,000V)

- **Mass spectrometry:** Orbitrap Exploris 480. MS1 settings: accumulation time = 250 ms (0-10 cells) or 100 ms (150-500 cells); resolution = 120,000; AGC = 100\ duration = 90 s (0-10 cells) or 60 s (150-500 cells) ; accumulation time = 500 ms (0-1 cell), 250 ms (10 cells), 100 ms (150 cells) or 50 ms (500 cells); resolution = 60,000 (0-10 cells) or 30,000 (150-500 cells); AGC = 5E3 (0-1 cells) or 1E4 (10-500 cells).
- **Data analysis:** MaxQuant (v1.6.7.0) and the search database is Swiss-Prot (July 2020).

### Data collection

All data were collected from the PRIDE repository (accession ID: PXD021882).

The sample annotations were collected from the methods section and from table S3 in the paper.

The PSM data were found in the evidence.txt file. The data were converted to a `QFeatures::QFeatures` object using the `scp::readSCP()` function.

The peptide data were found in the peptides.txt file. The column names holding the quantitative data were adapted to match the sample names in the `QFeatures::QFeatures` object. The data were then converted to a `SingleCellExperiment::SingleCellExperiment` object and then inserted in the `QFeatures::QFeatures` object. Links between the PSMs and the peptides were added

A similar procedure was applied to the protein data. The data were found in the proteinGroups.txt file. The column names were adapted, the data were converted to a `SingleCellExperiment::SingleCellExperiment` object and then inserted in the `QFeatures::QFeatures` object. Links between the peptides and the proteins were added

### Source

The PSM data can be downloaded from the PRIDE repository PXD021882 The source link is: <http://ftp.pride.ebi.ac.uk/pride/data/archive/2020/12/PXD021882/>

### References

Liang, Yiran, Hayden Acor, Michaela A. McCown, Andikan J. Nwosu, Hannah Boekweg, Nathaniel B. Axtell, Thy Truong, Yongzheng Cong, Samuel H. Payne, and Ryan T. Kelly. 2020. "Fully Automated Sample Processing and Analysis Workflow for Low-Input Proteome Profiling." *Analytical Chemistry*, December. ([link to article](#)).

### Examples

```
liang2020_hela()
```

---

petrosius2023\_AstralAML

*Petrosius et al. 2023 (bioRxiv): AML hierarchy on Astral.*

---

### Description

Single cell proteomics data from FACS sorted cells from the OCI-AML8227 model. The dataset contains leukemic stem cells (LSC; CD34+, CD38-), progenitor cells (CD34+, CD38+), CD38+ blasts (CD34-, CD38+) and CD38- blasts (CD34-, CD38-). It contains quantitative information at PSM, peptide and protein levels. Data was acquired using an Orbitrap Astral mass spectrometer. Direct DIA analysis was performed with Spectronaut version 17.

## Usage

petrosius2023\_AstralAML

## Format

A `QFeatures::QFeatures` object with 217 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- Assays 1-215: PSM data from the Spectronaut PEPQuant file with LFQ quantities from the FG.MS1Quantity column.
- peptides: Peptide data resulting from the PSM to peptide aggregation the 215 PSM assays. Resulting peptide assays were joined into a single assay.
- proteins: Protein data from the Spectronaut PGQuant file with LFQ quantities from the PG.Quantity column.

The `colData(petrosius2023_AstralAML())` contains cell type annotation, batch annotation and FACS data. The description of the `rowData` fields can be found in the [Spectronaut user manual](#).

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see *References*).

- **Cell isolation:** Cell sorting was done on a FACS Aria III or Aria II instrument, controlled by the DIVA software package and operated with a 100 microm nozzle. Cells were sorted at single-cell resolution, into a 384-well Eppendorf LoBind PCR plate containing 1 microL of lysis buffer.
- **Sample preparation** Single-cell protein lysates were digested overnight at 37°C with 2 ng of Trypsin supplied in 1 microL of digestion buffer. Digestion was stopped by the addition of 1 microL 1% (v/v) trifluoroacetic acid (TFA). All liquid dispensing was done using an I-DOT One instrument.
- **Liquid chromatography:** Chromatographic separation of peptides was conducted on a vanquish Neo UHPLC system connected to a 50 cm uPAC Neo Low-load and an EASY-spray. Autosampler and injection valves were configured to perform direct injections from a 384 well plate using a 25 uL injection loop on 11.8 min gradients.
- **Mass spectrometry:** Acquisition was conducted with an Orbitrap Astral mass spectrometer operated in positive mode with the FAIMSPro interface compensation voltage set to -45 V. MS1 scans were acquired with the Orbitrap at a resolution of 120,000 and a scan range of 400 to 900 m/z with normalized automatic gain control (AGC) target of 300 % and maximum injection time of 246 ms. Data independent acquisition of MS2 spectra was performed in the Astral using loop control set to 0.7 seconds per cycle with varying isolation window widths and injection times. Fragmentation of precursor ions was performed using higher energy collisional dissociation (HCD) using a normalized collision energy (NCE) of 25 %. AGC target was set to 800 %.
- **Raw data processing:** Raw files were processed using Spectronaut version 17. Direct DIA analysis was performed in pipeline mode. Pulsar searches were performed without fixed modifications. N-terminal acetylation and methionine oxidation were set as variable modifications. Quantification level was set to MS1 and the quantity type set to area under the curve.

## Data collection

The data were provided by the authors and is accessible at the [Dataverse](#). The dataset ('Astral AML single-cell data from Petrosius et al. 2023 preprint') contains the following files of interest:

- 20240201\_130747\_PEPQuant (Normal).tsv: the PSM level data
- 20240201\_130747\_PGQuant (Normal).tsv: the protein level data
- index\_map.csv: FACS data.
- msRuns\_overview.csv: Sample annotations.

We added the FACS data to the sample annotations in a single table. Both annotations and PSM features tables are then combined in a single `QFeatures::QFeatures` object using the `scp::readSCP()` function.

The peptide data were obtained by aggregation of the PSM data to the peptide level. All of the resulting peptides assays were joined into a single assays. Individual peptides assays were discarded.

The protein data were formatted from the 20240201\_130747\_PGQuant (Normal).tsv to a `SingleCellExperiment::SingleCellExperiment` object and the sample metadata were matched to the column names and stored in the `colData`. The object is then added to the `QFeatures::QFeatures` object and the rows of the peptide data are linked to the rows of the protein data based on the protein sequence information through an `AssayLink` object.

Note that the `QFeatures::QFeatures` object has not been further processed and has therefore not been normalized, log-transformed or batch-corrected.

## Source

The PSM data, protein data and sample annotations can be downloaded from the dataset 'Astral AML single-cell data from Petrosius et al. 2023 preprint' in the [Dataverse](#).

## References

Valdemaras Petrosius, Pedro Aragon-Fernandez, Tabiwang N. Arrey, Nil Üresin, Benjamin Furtwängler, Hamish Stewart, Eduard Denisov, Johannes Petzoldt, Amelia C. Peterson, Christian Hock, Eugen Damoc, Alexander Makarov, Vlad Zabrouskov, Bo T. Porse and Erwin M. Schoof. 2023. "Evaluating the capabilities of the Astral mass analyzer for single-cell proteomics." *bioRxiv*. <https://doi.org/10.1101/2023.06.06.543943>

## Examples

```
petrosius2023_AstralAML()
```

---

petrosius2023_mES	<i>Petrosius et al, 2023 (Nat. Comm.): Mouse embryonic stem cell (mESC) in different culture conditions</i>
-------------------	---

---

## Description

Profiling mouse embryonic stem cells across ground-state (m2i) and differentiation-permissive (m15) culture conditions. The data were acquired using orbitrap-based data-independent acquisition (DIA). The objective was to demonstrate the capability of their approach by profiling mouse embryonic stem cell culture conditions, showcasing heterogeneity in global proteomes, and highlighting differences in the expression of key metabolic enzymes in distinct cell subclusters.

## Usage

petrosius2023\_mES

## Format

A `QFeatures::QFeatures` object with 605 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- Assay 1-603: PSM data acquired with an orbitrap-based data-independent acquisition (DIA) protocol, hence those assays contain single column that contains the quantitative information.
- peptides: peptide data containing quantitative data for 9884 peptides and 603 single-cells.
- proteins: protein data containing quantitative data for 4270 proteins and 603 single-cells.

Sample annotation is stored in `colData(petrosius2023_mES())`.

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Sample isolation:** Cell sorting was done on a Sony MA900 cell sorter using a 130 microm sorting chip. Cells were sorted at single-cell resolution, into a 384-well Eppendorf LoBind PCR plate (Eppendorf AG) containing 1 microL of lysis buffer.
- **Sample preparation:** Single-cell protein lysates were digested with 2 ng of Trypsin supplied in 1 microL of digestion buffer which was carried out overnight at 37 °C, and subsequently acidified by the addition of 1 microL 1% (v/v) trifluoroacetic acid (TFA). All liquid dispensing was done using an I-DOT One instrument.
- **Liquid chromatography:** For the HRMS1-DIA experiments and the DIA isolation window survey, Evosep One liquid chromatography was used. The standard 31-minute or 58-minute pre-defined Whisper gradients were used with a flow rate of 100 nl/min for peptide elution.
- **Mass spectrometry:** The mass spectrometer was operated in positive mode with the FAIM-SPro interface compensation voltage set to -45 V. MS1 scans were carried out at 120,000 resolution with an automatic gain control (AGC) of 300% and maximum injection time set to auto. For the DIA isolation window survey, a scan range of 500–900 was used, and 400–1000 rest of the experiments. Higher energy collisional dissociation (HCD) was used for precursor fragmentation with a normalized collision energy (NCE) of 33% and the MS2 scan AGC target was set to 1000%.
- **Raw data processing:** The mESC raw data files were processed with Spectronaut 17.

## Data collection

The data were provided by the Author and is accessible at the [Dataverse](#) The folder ('20240205\_111248\_mESC\_SNEcom m2i/') contains the following files of interest:

- 20240205\_111251\_PEPQuant (Normal).tsv: the PSM level data
- 20240205\_111251\_Peptide Quant (Normal).tsv: the peptide level data
- 20240205\_111251\_PGQuant (Normal).tsv: the protein level data

The metadata were downloaded from the [Zenodo repository](#).

- sample\_facs.csv: the metadata

We formatted the quantification table so that columns match with the metadata. Then, both tables are then combined in a single `QFeatures::QFeatures` object using the `scp::readSCP()` function.

The peptide data were formatted to a `SingleCellExperiment::SingleCellExperiment` object and the sample metadata were matched to the column names and stored in the `colData`. The object is then added to the `QFeatures::QFeatures` object and the rows of the PSM data are linked to the rows of the peptide data based on the peptide sequence information through an `AssayLink` object.

The protein data were formatted to a `SingleCellExperiment::SingleCellExperiment` object and the sample metadata were matched to the column names and stored in the `colData`. The object is then added to the `QFeatures::QFeatures` object and the rows of the peptide data are linked to the rows of the protein data based on the protein sequence information through an `AssayLink` object.

## Source

The peptide and protein data can be downloaded from the [Dataverse](#). The raw data and the quantification data can also be found in the MassIVE repository `MSV000092429`: `ftp://MSV000092429@massive.ucsd.edu/`.

## References

**Source article:** Petrosius, V., Aragon-Fernandez, P., Üresin, N. et al. "Exploration of cell state heterogeneity using single-cell proteomics through sensitivity-tailored data-independent acquisition." *Nat Commun* 14, 5910 (2023). ([link to article](#)).

## Examples

```
petrosius2023_mES()
```

---

schoof2021

*Schoof et al. 2021 (Nat. Comm.): acute myeloid leukemia differentiation*

---

## Description

Single-cell proteomics data from OCI-AML8227 cell culture to reconstruct the cellular hierarchy. The data were acquired using TMTpro multiplexing. The samples contain either no cells, single cells, 10 cells (reference channel) 200 cells (booster channel) or are simply empty wells. Single cells are expected to be one of progenitor cells (PROG), leukaemia stem cells (LSC), CD38- blast cells (BLAST CD38-) or CD38+ blast cells (BLAST CD38+). Booster are either a known 1:1:1 mix of cells (PROG, LSC and BLAST) or are isolated directly from the bulk sample. Samples were isolated and annotated using flow cytometry.

## Usage

```
schoof2021
```

## Format

A `QFeatures::QFeatures` object with 194 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `F*`: 192 assays containing PSM quantification data for 16 TMT channels. The quantification data contain signal to noise ratios as computed by Proteome Discoverer.
- `proteins`: quantitative data for 2898 protein groups in 3072 samples (all runs combined). The quantification data contain signal to noise ratios as computed by Proteome Discoverer.
- `logNormProteins`: quantitative data for 2723 protein groups in 2025 single-cell samples. This assay is the protein datasets that was processed by the authors. Dimension reduction and clustering data are also available in the `reducedDims` and `colData` slots, respectively

Sample annotation is stored in `colData(schoof2021())`. The cell type annotation is stored in the `Population` column. The flow cytometry data is also available: `FSC-A`, `FSC-H`, `FSC-W`, `SSC-A`, `SSC-H`, `SSC-W`, `APC-Cy7-A` (= CD34) and `PE-A` (= CD38).

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Sample isolation**: cultured AML 8227 cells were stained with anti-CD34 and anti-CD38. The sorting was performed by FACS Aria instrument and deposited in 384 well plates.
- **Sample preparation**: cells are lysed using freeze-boil and sonication in a lysis buffer (TFE) that also includes reduction and alkylation reagents (TCEP and CAA), followed by trypsin (protein) and benzonase (DNA) digestion, TMT-16 labeling and quenching, desalting using SOLA $\mu$  C18 plate, peptide concentration, pooling and peptide concentration again. The booster channel contains 200 cell equivalents.
- **Liquid chromatography**: peptides are separated using a C18 reverse-phase column (50cm x 75  $\mu$ m i.d., Thermo EasySpray) combined to a Thermo EasyLC 1200 for 160 minute gradient with a flowrate of 100nl/min.
- **Mass spectrometry**: FAIMSPro interface is used. MS1 setup: resolution 60.000, AGC target of 300%, accumulation of 50ms. MS2 setup: resolution 45.000, AGC target of 150, 300 or 500%, accumulation of 150, 300, 500, or 1000ms.
- **Raw data processing**: Proteome Discoverer 2.4 + Sequest spectral search engine and validation with Percolator

## Data collection

All data were collected from the PRIDE repository (accession ID: PXD020586). The data and metadata were extracted from the `SCEptre_FINAL.zip` file.

We performed extensive data wrangling to combine all the metadata available from different files into a single table available using `colData(schoof2021)`.

The PSM data were found in the `bulk_PSMs.txt` file. Contaminants were defined based on the protein accessions listed in `contaminant.txt`. The data were converted to a `QFeatures::QFeatures` object using the `scp::readSCP()` function.

The protein data were found in the `bulk_Proteins.txt` file. Contaminants were defined based on the protein accessions listed in `contaminant.txt`. The column names holding the quantitative data were adapted to match the sample names in the `QFeatures::QFeatures` object. Unnecessary feature annotations (such as in which assay a protein is found) were removed. Feature

names were created following the procedure in SCEPTRE: features names are the protein symbol (or accession if missing) and if duplicated symbols are present (protein isoforms), they are made unique by appending the protein accession. Contaminants were defined based on the protein accessions listed in `contaminant.txt`. The data were then converted to a `SingleCellExperiment::SingleCellExperiment` object and inserted in the `QFeatures::QFeatures` object.

The log-normalized protein data were found in the `bulk.h5ad` file. This dataset was generated by the authors by running the notebook called `bulk.ipynb`. The `bulk.h5ad` was loaded as an `AnnData` object using the `scanpy` Python module. The object was then converted to a `SingleCellExperiment` object using the `zellkonverter` package. The column names holding the quantitative data were adapted to match the sample names in the `QFeatures::QFeatures` object. The data were then inserted in the `QFeatures::QFeatures` object.

The script to reproduce the `QFeatures` object is available at `system.file("scripts", "make-data_schoof2021.R", package = "scpdata")`

### Source

The PSM and protein data can be downloaded from the PRIDE repository PXD020586 The source link is: <https://www.ebi.ac.uk/pride/archive/projects/PXD020586>

### References

Schoof, Erwin M., Benjamin Furtwängler, Nil Üresin, Nicolas Rapin, Simonas Savickas, Coline Gentil, Eric Lechman, Ulrich auf Dem Keller, John E. Dick, and Bo T. Porse. 2021. "Quantitative Single-Cell Proteomics as a Tool to Characterize Cellular Hierarchies." *Nature Communications* 12 (1): 745679. ([link to article](#)).

### Examples

```
schoof2021()
```

---

scpdata

*Single-Cell Proteomics Data Package*

---

### Description

The `scpdata` package distributes mass spectrometry-based single-cell proteomics datasets. The datasets were collected from published work and formatted to a standardized data framework. The `scpdata` framework stores the expression data for different MS levels (identified spectrum, peptide, or protein) in separate assays. Each assay is an object of class `SingleCellExperiment::SingleCellExperiment` that allows easy integration with state-of-the-art single-cell analysis tools. All assays are contained in a single object of class `QFeatures::QFeatures`. An overview of the data structure is shown provided in the `scpdata` package.

The `scpdata()` function returns a summary table with all currently available datasets in the package. More information about the data content and the data collection can be found in the corresponding manual pages.

### Usage

```
scpdata()
```

**Value**

A DataFrame table containing a summary of the available datasets.

**Author(s)**

Christophe Vanderaa

**See Also**

More information about the data manipulation can be found in the `scp` package.

**Examples**

```
## List available datasets and their metadata
scpdata()

## Load data using the ExperimentHub interface
hub <- ExperimentHub()

## Not run:
## Download the data set of interest using ExperimentHub indexing
hub[["EH7711"]]
## Download the same data set using the build-in function
leduc2022()

## End(Not run)
```

---

specht2019v2

*Specht et al. 2019 - SCoPE2 (biorRxiv): macrophages vs monocytes (version 2)*

---

**Description**

Single cell proteomics data acquired by the Slavov Lab. This is the version 2 of the data released in December 2019. It contains quantitative information of macrophages and monocytes at PSM, peptide and protein level.

**Usage**

```
specht2019v2
```

**Format**

A `QFeatures::QFeatures` object with 179 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- Assay 1-63: PSM data for SCoPE2 sets acquired with a TMT-11plex protocol, hence those assays contain 11 columns. Columns hold quantitative information from single-cell channels, carrier channels, reference channels, empty (blank) channels and unused channels.
- Assay 64-177: PSM data for SCoPE2 sets acquired with a TMT-16plex protocol, hence those assays contain 16 columns. Columns hold quantitative information from single-cell channels, carrier channels, reference channels, empty (blank) channels and unused channels.

- peptides: peptide data containing quantitative data for 9208 peptides and 1018 single-cells.
- proteins: protein data containing quantitative data for 2772 proteins and 1018 single-cells.

The `colData(specht2019v2())` contains cell type annotation and batch annotation that are common to all assays. The description of the `rowData` fields for the PSM data can be found in the [MaxQuant documentation](#).

### Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation:** flow cytometry (BD FACSAria I).
- **Sample preparation** performed using the SCoPE2 protocol. mPOP cell lysis + trypsin digestion + TMT-11plex or 16plex labelling and pooling.
- **Separation:** online nLC (DionexUltiMate 3000 UHPLC with a 25cm x 75um IonOpticksAurora Series UHPLC column; 200nL/min).
- **Ionization:** ESI (2,200V).
- **Mass spectrometry:** Thermo Scientific Q-Exactive (MS1 resolution = 70,000; MS1 accumulation time = 300ms; MS2 resolution = 70,000).
- **Data analysis:** DART-ID + MaxQuant (1.6.2.3).

### Data collection

The PSM data were collected from a shared Google Drive folder that is accessible from the SlavovLab website (see Source section). The folder contains the following files of interest:

- `ev_updated.txt`: the MaxQuant/DART-ID output file
- `annotation_fp60-97.csv`: sample annotation
- `batch_fp60-97.csv`: batch annotation

We combined the sample annotation and the batch annotation in a single table. We also formatted the quantification table so that columns match with those of the annotation and filter only for single-cell runs. Both table are then combined in a single `QFeatures::QFeatures` object using the `scp::readSCP()` function.

The peptide data were taken from the Slavov lab directly (`Peptides-raw.csv`). It is provided as a spreadsheet. The data were formatted to a `SingleCellExperiment::SingleCellExperiment` object and the sample metadata were matched to the column names (mapping is retrieved after running the SCoPE2 R script) and stored in the `colData`. The object is then added to the `QFeatures::QFeatures` object (containing the PSM assays) and the rows of the peptide data are linked to the rows of the PSM data based on the peptide sequence information through an `AssayLink` object.

The protein data (`Proteins-processed.csv`) is formatted similarly to the peptide data, and the rows of the proteins were mapped onto the rows of the peptide data based on the protein sequence information.

### Source

The data were downloaded from the [Slavov Lab](#) website via a shared Google Drive [folder](#). The raw data and the quantification data can also be found in the massIVE repository `MSV000083945`: `ftp://massive.ucsd.edu/MSV000083945`.

## References

Specht, Harrison, Edward Emmott, Aleksandra A. Petelski, R. Gray Huffman, David H. Perlman, Marco Serra, Peter Kharchenko, Antonius Koller, and Nikolai Slavov. 2019. "Single-Cell Mass-Spectrometry Quantifies the Emergence of Macrophage Heterogeneity." bioRxiv. ([link to article](#)).

## Examples

specht2019v2()

---

specht2019v3	<i>Specht et al. 2019 - SCoPE2 (biorRxiv): macrophages vs monocytes (version 3)</i>
--------------	---

---

## Description

Single cell proteomics data acquired by the Slavov Lab. This is the version 3 of the data released in October 2020. It contains quantitative information of macrophages and monocytes at PSM, peptide and protein level.

## Usage

specht2019v3

## Format

A `QFeatures::QFeatures` object with 179 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- Assay 1-63: PSM data for SCoPE2 sets acquired with a TMT-11plex protocol, hence those assays contain 11 columns. Columns hold quantitative information from single-cell channels, carrier channels, reference channels, empty (blank) channels and unused channels.
- Assay 64-177: PSM data for SCoPE2 sets acquired with a TMT-16plex protocol, hence those assays contain 16 columns. Columns hold quantitative information from single-cell channels, carrier channels, reference channels, empty (blank) channels and unused channels.
- peptides: peptide data containing quantitative data for 9208 peptides and 1018 single-cells.
- proteins: protein data containing quantitative data for 2772 proteins and 1018 single-cells.

The `colData(specht2019v2())` contains cell type annotation and batch annotation that are common to all assays. The description of the `rowData` fields for the PSM data can be found in the [MaxQuant documentation](#).

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation:** flow cytometry (BD FACSAria I).
- **Sample preparation** performed using the SCoPE2 protocol. mPOP cell lysis + trypsin digestion + TMT-11plex or 16plex labeling and pooling.

- **Separation:** online nLC (DionexUltiMate 3000 UHPLC with a 25cm x 75um IonOptick-sAurora Series UHPLC column; 200nL/min).
- **Ionization:** ESI (2,200V).
- **Mass spectrometry:** Thermo Scientific Q-Exactive (MS1 resolution = 70,000; MS2 accumulation time = 300ms; MS2 resolution = 70,000).
- **Data analysis:** DART-ID + MaxQuant (1.6.2.3).

## Data collection

The PSM data were collected from a shared Google Drive folder that is accessible from the SlavovLab website (see Source section). The folder contains the following files of interest:

- `ev_updated_v2.txt`: the MaxQuant/DART-ID output file
- `annotation_fp60-97.csv`: sample annotation
- `batch_fp60-97.csv`: batch annotation

We combined the sample annotation and the batch annotation in a single table. We also formatted the quantification table so that columns match with those of the annotation and filter only for single-cell runs. Both table are then combined in a single `QFeatures::QFeatures` object using the `scp::readSCP()` function.

The peptide data were taken from the Slavov lab directly (`Peptides-raw.csv`). It is provided as a spreadsheet. The data were formatted to a `SingleCellExperiment::SingleCellExperiment` object and the sample metadata were matched to the column names (mapping is retrieved after running the `SCoPE2 R` script) and stored in the `colData`. The object is then added to the `QFeatures::QFeatures` object (containing the PSM assays) and the rows of the peptide data are linked to the rows of the PSM data based on the peptide sequence information through an `AssayLink` object.

The protein data (`Proteins-processed.csv`) is formatted similarly to the peptide data, and the rows of the proteins were mapped onto the rows of the peptide data based on the protein sequence information.

## Note

Since version 2, a serious bug in the data were corrected for TMT channels 12 to 16. Many more cells are therefore contained in the data. Version 2 is maintained for backward compatibility. Although the final version of the article was published in 2021, we have kept `specht2019v3` as the data set name for consistency with the previous data version `specht2019v2`.

## Source

The data were downloaded from the [Slavov Lab](#) website via a shared Google Drive [folder](#). The raw data and the quantification data can also be found in the massIVE repository `MSV000083945`: <ftp://massive.ucsd.edu/MSV000083945>.

## References

Specht, Harrison, Edward Emmott, Aleksandra A. Petelski, R. Gray Huffman, David H. Perlman, Marco Serra, Peter Kharchenko, Antonius Koller, and Nikolai Slavov. 2021. "Single-Cell Proteomic and Transcriptomic Analysis of Macrophage Heterogeneity Using SCoPE2." *Genome Biology* 22 (1): 50. ([link to article](#)).

## Examples

```
specht2019v3()
```

---

```
williams2020_lfq
```

```
Williams et al. 2020 (Anal. Chem.): MCF10A cell line
```

---

## Description

Single-cell label free proteomics data from a MCF10A cell line culture. The data were acquired using a label-free quantification protocol based on the nanoPOTS technology. The objective was to test 2 elution gradients for single-cell applications and to demonstrate successful use of the new nanoPOTS autosampler presented in the article. The samples contain either no cells, single cells, 3 cells, 10 cells 50 cells.

## Usage

```
williams2020_lfq
```

## Format

A `QFeatures::QFeatures` object with 9 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `peptides_[30 or 60]min_[intensity or LFQ]`: 3 assays containing peptide intensities or LFQ normalized quantifications (see References) for either a 30min or a 60 min gradient.
- `proteins_[30 or 60]min_[intensity or iBAQ or LFQ]`: 6 assays containing protein intensities, iBAQ normalized or LFQ normalized quantifications (see References) for either a 30min or a 60 min gradient.

Sample annotation is stored in `colData(williams2020_lfq())`.

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Sample isolation:** cultured MCF10A cells were isolated using flow-cytometry based cell sorting and deposit on nanoPOTS microwells
- **Sample preparation:** cells are lysed using using a DDM+DTT lysis buffer. Alkylation was then performed using an IAA solution. Proteins are digested with Lys-C and trypsin followed by acidification with FA. Sample droplets are then dried until LC-MS/MS analysis.
- **Liquid chromatography:** peptides are loaded using the new autosampler described in the paper. Samples are loaded using a a homemade miniature syringe pump. The samples are then desalted and concentrated through a SPE column (4cm x 100µm i.d. packed with 5µm C18) with microflow LC pump. The peptides are then eluted from a long LC column (60cm x 50 µm i.d. packed with 3µm C18) coupled to a nanoflow LC pump at 150nL/mL with either a 30 min or a 60 min gradient.

- **Mass spectrometry:** MS/MS was performed on an Orbitrap Fusion Lumos Tribrid MS coupled to a 2kV ESI. MS1 setup: Orbitrap analyzer at resolution 120.000, AGC target of 1E6, accumulation of 246ms. MS2 setup: ion trap with CID at resolution 60.000, AGC target of 2E4, accumulation of 120ms (50 cells) or 250ms (0-10 cells).
- **Raw data processing:** preprocessing using Maxquant v1.6.2.10 that use Andromeda search engine (with UniProtKB 2016-21-29), MBR and LFQ normalization were enabled.

### Data collection

All data were collected from the MASSIVE repository (accession ID: MSV000085230).

The peptide and protein data were extracted from the Peptides\_[...].txt or ProteinGroups\_[...].txt files, respectively, in the MCF10A\_LC\_[30 or 60]minutes folders.

The tables were duplicated so that peptide intensities, peptide LFQ, protein intensities, protein LFQ and protein intensities are contained in separate tables. Tables are then converted to `SingleCellExperiment::SingleCellExperiment` objects. Sample annotations were inferred from the sample names and from the paper. All data is combined in a `QFeatures::QFeatures` object. `QFeatures::AssayLinks` were stored between peptide assays and their corresponding proteins assays based on the leading razor protein (hence only unique peptides are linked to proteins).

The script to reproduce the QFeatures object is available at `system.file("scripts", "make-data_williams2020_lfq", package = "scpdata")`

### Suggestion

See `QFeatures::joinAssays` if you want to join the 30min and 60min assays in a single assay for an integrated analysis.

### Source

The PSM and protein data can be downloaded from the MASSIVE repository MSV000085230.

### References

**Source article:** Williams, Sarah M., Andrey V. Liyu, Chia-Feng Tsai, Ronald J. Moore, Daniel J. Orton, William B. Chrisler, Matthew J. Gaffrey, et al. 2020. "Automated Coupling of Nanodroplet Sample Preparation with Liquid Chromatography-Mass Spectrometry for High-Throughput Single-Cell Proteomics." *Analytical Chemistry* 92 (15): 10588–96. ([link to article](#)).

**LFQ normalization:** Cox, Jürgen, Marco Y. Hein, Christian A. Lubner, Igor Paron, Nagarjuna Nagaraj, and Matthias Mann. 2014. "Accurate Proteome-Wide Label-Free Quantification by Delayed Normalization and Maximal Peptide Ratio Extraction, Termed MaxLFQ." *Molecular & Cellular Proteomics: MCP* 13 (9): 2513–26. ([link to article](#)).

**iBAQ normalization:** Schwanhäusser, Björn, Dorothea Busse, Na Li, Gunnar Dittmar, Johannes Schuchhardt, Jana Wolf, Wei Chen, and Matthias Selbach. 2011. "Global Quantification of Mammalian Gene Expression Control." *Nature* 473 (7347): 337–42. ([link to article](#)).

### Examples

```
williams2020_lfq()
```

---

williams2020\_tmt

*Williams et al. 2020 (Anal. Chem.): 3 AML cell line*

---

## Description

Single-cell label data from three acute myeloid leukemia cell line culture (MOLM-14, K562, CMK). The data were acquired using a TMT-based quantification protocol and the nanoPOTS technology. The objective was to demonstrate successful use of the new nanoPOTS autosampler presented in the source article. The samples contain either carrier (10 ng), reference (0.2ng), empty or single-cell samples..

## Usage

williams2020\_tmt

## Format

A `QFeatures::QFeatures` object with 4 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `peptides_[intensity or corrected]`: 2 assays containing peptide reporter ion intensities or corrected reporter ion intensities as computed by MaxQuant.
- `proteins_[intensity or corrected]`: 2 assays containing protein reporter ion intensities or corrected reporter ion intensities as computed by MaxQuant.

Sample annotation is stored in `colData(williams2020_tmt())`.

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Sample isolation:** cultured MOLM-14, K562 or CMK cells were isolated using flow-cytometry based cell sorting and deposit on nanoPOTS microwells
- **Sample preparation:** cells are lysed using using a DDM lysis buffer. Proteins are digested with trypsin followed by TMT labelling and quanching with HA. The samples are then acidified with FA, pooled in a single samples (adding carrier and reference peptide mixtures), and dried until LC-MS/MS analysis.
- **Liquid chromatography:** peptides are loaded using the new autosampler described in the paper. Samples are loaded using a a homemade miniature syringe pump. The samples are then desalted and concentrated through a SPE column (4cm x 100µm i.d. packed with 5µm C18) with microflow LC pump. The peptides are then eluted from a long LC column (60cm x 50 µm i.d. packed with 3µm C18) coupled to a nanoflox LC pump at 150nL/mL (elution time is not expliceted).
- **Mass spectrometry:** MS/MS was performed on an Orbitrap Fusion Lumos Tribrid MS coupled to a 2kV ESI. MS1 setup: Orbitrap analyzer at resolution 120.000, AGC target of 1E6, accumulation of 246ms. MS2 setup: Orbitrap with HCD at resolution 120.000, AGC target of 1E6, accumulation of 246ms.
- **Raw data processing:** preprocessing using Maxquant v1.6.2.10 that use Andromeda search engine (with UniProtKB 2016-21-29).

## Data collection

All data were collected from the MASSIVE repository (accession ID: MSV000085230).

The peptide and protein data were extracted from the `Peptides_AML_SingleCell.txt` or `ProteinGroups_AML_SingleCell.txt` files, respectively, in the `AML_SingleCell` folders.

The tables were duplicated so that intensities and corrected intensities are contained in separate tables. Tables are then converted to `SingleCellExperiment::SingleCellExperiment` objects. Sample annotations were inferred from the sample names, from table S2 and from the Experimental Section of the paper. All data is combined in a `QFeatures::QFeatures` object. `QFeatures::AssayLinks` were stored between peptide assays and their corresponding proteins assays based on the leading razor protein (hence only unique peptides are linked to proteins).

The script to reproduce the `QFeatures` object is available at `system.file("scripts", "make-data_williams2020_tmt", package = "scpdata")`

## Source

The PSM and protein data can be downloaded from the MASSIVE repository MSV000085230.

## References

**Source article:** Williams, Sarah M., Andrey V. Liyu, Chia-Feng Tsai, Ronald J. Moore, Daniel J. Orton, William B. Chrisler, Matthew J. Gaffrey, et al. 2020. "Automated Coupling of Nanodroplet Sample Preparation with Liquid Chromatography-Mass Spectrometry for High-Throughput Single-Cell Proteomics." *Analytical Chemistry* 92 (15): 10588–96. ([link to article](#)).

## Examples

```
williams2020_tmt()
```

---

woo2022\_lung

*Woo et al. 2022 (Cell Syst.): 26 primary human lung cells*

---

## Description

Single-cell proteomics data from dissociated primary human lung cells. The data were acquired using the TIFF (transfer identification based on FAIMS filtering) acquisition method. The data contain 26 single cells.

## Usage

```
woo2022_lung
```

## Format

A `QFeatures::QFeatures` object with 5 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `peptides_[intensity or LFQ]`: 2 assays containing peptide quantities or normalized quantities using the maxLFQ method as computed by `MaxQuant`.

- `proteins_[intensity or iBAQ or LFQ]`: 3 assays containing protein quantities or normalized proteins using the iBAQ or maxLFQ methods as computed by MaxQuant.

Sample annotation is stored in `colData(woo_lung())`.

### Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Sample isolation:** primary human lung cells were dissociated following the protocol in Bandyopadhyay et al., 2018. The cells were sorted using the Influx II cell sorter and deposited on a nanoPOTS chip.
- **Sample preparation:** cells are lysed using using a DDM+DTT lysis and reduction buffer. The proteins are alkylated with IAA and digested with LysC and trypsin. Samples are then acidified with FA, vacuum dried and stored in freezer until data acquisition.
- **Liquid chromatography:** peptides are loaded using an in-house autosampler (Williams et al. 2020). The samples are concentrated through a SPE column (4cm x 100µm i.d. packed with 5µm C18) with microflow LC pump. The peptides are then eluted from an LC column (25cm x 50 µm i.d. packed with 1.7µm C18) from a 60 min gradient (100nL/min).
- **Mass spectrometry:** MS/MS was performed on an Orbitrap Fusion Lumos Tribrid MS with FAIMSpro coupled to a 2.4 kV ESI. FAIMS setup: 4-CV method (-45, -55, -65, -75 V). MS1 setup: resolution = 120,000, range = 350-1500 m/z, AGC target of 1E6, accumulation of 254ms. MS2 setup: 30% HCD, resolution AGC 2E4, accumulation of 254ms.
- **Raw data processing:** preprocessing using Maxquant v1.6.2.10 that use Andromeda search engine (with UniProtKB 2016-21-29). MBR was enabled.

### Data collection

All data were collected from the MASSIVE repository (accession ID: MSV000085937).

The peptide and protein data were extracted from the `peptides_nondepleted_Lung_scProteomics.txt` or `proteinGroups_nondepleted_Lung_scProteomics.txt` files, respectively, in the `NonDepleted_Lung_SingleCell` folders.

The tables were split so that intensities, maxLFQ, and iBAQ data are contained in separate tables. Tables are then converted to `SingleCellExperiment::SingleCellExperiment` objects. Sample annotations were inferred from the sample names. All data is combined in a `QFeatures::QFeatures` object. `QFeatures::AssayLinks` were stored between peptide assays and their corresponding proteins assays based on the leading razor protein (hence only unique peptides are linked to proteins).

The script to reproduce the `QFeatures` object is available at `system.file("scripts", "make-data_woo2022_lung.R")` package = "scpdata")

### Source

The peptide and protein data can be downloaded from the MASSIVE repository MSV000085937

### References

**Source article:** Woo, Jongmin, Jeremy C. Clair, Sarah M. Williams, Song Feng, Chia-Feng Tsai, Ronald J. Moore, William B. Chrisler, et al. 2022. "Three-Dimensional Feature Matching Improves Coverage for Single-Cell Proteomics Based on Ion Mobility Filtering." *Cell Systems* 13 (5): 426–34.e4. ([link to article](#)).

## Examples

```
woo2022_lung()
```

---

woo2022_macrophage	<i>Woo et al. 2022 (Cell Syst.): LPS-treated macrophages</i>
--------------------	--

---

## Description

Single-cell data from macrophages subjected to 3 LPS treatments. The data were acquired using the TIFF (transfer identification based on FAIMS filtering) acquisition method. The data contain 155 single cells: 54 control cells (no treatment), 52 cells treated with LPS during 24h and 49 cells treated with LPS during 49h.

## Usage

```
woo2022_macrophage
```

## Format

A `QFeatures::QFeatures` object with 5 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `peptides_[intensity or LFQ]`: 2 assays containing peptide quantities or normalized quantities using the maxLFQ method as computed by MaxQuant.
- `proteins_[intensity or iBAQ or LFQ]`: 3 assays containing protein quantities or normalized proteins using the iBAQ or maxLFQ methods as computed by MaxQuant.

Sample annotation is stored in `colData(woo_macrophage())`.

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Sample isolation:** cultured RAW 264.7 cells treated or not with 100 ng/ul LPS. The cells were sorted using the Influx II cell sorter and deposited on a nanoPOTS chip.
- **Sample preparation:** cells are lysed using using a DDM+DTT lysis and reduction buffer. The proteins are alkylated with IAA and digested with LysC and trypsin. Samples are then acidified with FA, vacuum dried and stored in freezer until data acquisition.
- **Liquid chromatography:** peptides are loaded using an in-house autosampler (Williams et al. 2020). The samples are concentrated through a SPE column (4cm x 100µm i.d. packed with 5µm C18) with microflow LC pump. The peptides are then eluted from an LC column (25cm x 50 µm i.d. packed with 1.7µm C18) from a 60 min gradient (100nL/min).
- **Mass spectrometry:** MS/MS was performed on an Orbitrap Fusion Lumos Tribrid MS with FAIMSpro coupled to a 2.4 kV ESI. FAIMS setup: 4-CV method (-45, -55, -65, -75 V). MS1 setup: resolution = 120.000, range = 350-1500 m/z, AGC target of 1E6, accumulation of 254ms. MS2 setup: 30% HCD, resolution AGC 2E4, accumulation of 254ms.
- **Raw data processing:** preprocessing using Maxquant v1.6.2.10 that use Andromeda search engine (with UniProtKB 2016-21-29). MBR was enabled.

## Data collection

All data were collected from the MASSIVE repository (accession ID: MSV000085937).

The peptide and protein data were extracted from the `peptides_RAW_LPS_scProteomics.txt` or `proteinGroups_RAW_LPS_scProteomics.txt` files, respectively, in the `RAW_LPS_SingleCellProteomics` folders.

The tables were split so that intensities, maxLFQ, and iBAQ data are contained in separate tables. Tables are then converted to `SingleCellExperiment::SingleCellExperiment` objects. Sample annotations were inferred from the sample names. All data is combined in a `QFeatures::QFeatures` object. `QFeatures::AssayLinks` were stored between peptide assays and their corresponding proteins assays based on the leading razor protein (hence only unique peptides are linked to proteins).

The script to reproduce the `QFeatures` object is available at `system.file("scripts", "make-data_woo2022_macrophage", package = "scpdata")`

## Source

The peptide and protein data can be downloaded from the MASSIVE repository MSV000085937

## References

**Source article:** Woo, Jongmin, Jeremy C. Clair, Sarah M. Williams, Song Feng, Chia-Feng Tsai, Ronald J. Moore, William B. Chrisler, et al. 2022. "Three-Dimensional Feature Matching Improves Coverage for Single-Cell Proteomics Based on Ion Mobility Filtering." *Cell Systems* 13 (5): 426–34.e4. ([link to article](#)).

## Examples

```
woo2022_macrophage()
```

---

wu2026

*Wu et al., 2026 (Nature Biotechnology): single-cell proteomic landscape of the developing human brain*

---

## Description

Label-free single-cell proteomics data from prenatal human brain tissue (gestational weeks 13, 15 and 19), generated to characterize cell-type proteomes and developmental trajectories at single-cell resolution.

## Usage

```
wu2026
```

## Format

A `QFeatures::QFeatures` object with 2313 sets, each set being a `SingleCellExperiment::SingleCellExperiment` object:

- Set 1-2310: log-transformed precursor-level single-cell data (one column per set) from DIA-NN reports.

- **peptides:** peptide-level data containing quantitative data for 17755 peptides in 2310 single cells.
- **imported\_proteins:** protein-level table imported from the authors' protein report, containing 4401 proteins in 2310 single cells.
- **proteins:** protein-level data containing quantitative data for 5107 proteins in 2310 single cells.

Sample annotation is stored in `colData(wu2026())`.

### Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Sample isolation:** fresh prenatal human brain tissue was dissociated and single cells were sorted by FACS (Sony SH800) into 384-well low-bind plates.
- **Sample preparation:** cells were lysed in low-volume TEAB-based buffer, digested overnight with trypsin, acidified with trifluoroacetic acid, and desalted on EvoTip Pure.
- **Liquid chromatography:** peptides were separated on an IonOpticks Aurora Elite XT C18 nano-LC column (15 cm x 75 micrometer) on an Evosep One using a 31-minute gradient.
- **Mass spectrometry:** Orbitrap Eclipse Tribrid with FAIMS Pro was used in DIA mode.
- **Raw data processing:** DIA-NN against the UniProt human reference proteome.

### Data collection

The `QFeatures` object was built from the DIA-NN peptide and protein report tables and a sample annotation table. PSM-level data were imported with `scp::readSCP()`, log-transformed, aggregated to peptide-level features, and joined across single-cell runs. The protein report table was imported as a `SingleCellExperiment::SingleCellExperiment`, added as a set, and linked to peptide-level data.

The script to reproduce the `QFeatures` object is available at `inst/scripts/make-data_wu2026.R`.

### Source

The data were downloaded from the PRIDE repository with accession ID PXD071075: <https://www.ebi.ac.uk/pride/archive>

### References

Wu, T., Jiang, L., Mukhtar, T., Wang, L., Jian, R., Wang, C., Trinh, T., Kriegstein, A. R., Snyder, M., and Li, J. 2026. "Single-cell proteomic landscape of the developing human brain." *Nature Biotechnology*. [Link to article](#)

### Examples

`wu2026()`

---

zhu2018MCP

*Zhu et al. 2018 (Mol. Cel. Prot.): rat brain laser dissections*

---

## Description

Near single-cell proteomics data of laser captured micro-dissection samples. The samples are 24 brain sections from rat pups (day 17). The slices are 12 um thick squares of either 50, 100, or 200 um width. 5 samples were dissected from the corpus callum (CC), 4 samples were dissected from the corpus collosum (CP), 13 samples were extracted from the cerebral cortex (CTX), and 2 samples are labeled as (Mix).

## Usage

zhu2018MCP

## Format

A `QFeatures::QFeatures` object with 4 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `peptides`: quantitative information for 13,055 peptides from 24 samples
- `proteins_intensity`: protein intensities for 2,257 proteins from 24 samples
- `proteins_LFQ`: LFQ intensities for 2,257 proteins from 24 samples
- `proteins_iBAQ`: iBAQ values for 2,257 proteins from 24 samples

Sample annotation is stored in `colData(zhu2018MCP())`.

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the original article (see References).

- **Cell isolation**: brain patches were collected using laser-capture microdissection (PALM MicroBeam) on flash frozen rat (*Rattus norvegicus*) brain tissues. Note that the samples were stained with H&E before dissection for histological analysis. DMSO is used as sample collection solution
- **Sample preparation** performed using the nanoPOTs device: DMSO evaporation + protein extraction (DMM + DTT) + alkylation (IAA)
  - Lys-C digestion + trypsin digestion.
- **Separation**: nanoLC (Dionex UltiMate with an in-house packed 60cm x 30um LC columns; 50nL/min)
- **Ionization**: ESI (2,000V)
- **Mass spectrometry**: Thermo Fisher Orbitrap Fusion Lumos Tribrid (MS1 accumulation time = 246ms; MS1 resolution = 120,000; MS1 AGC = 3E6). The MS/MS settings depend on the sample size, excepted for the AGC = 1E5. 50um (time = 502ms; resolution = 240,000), 100um (time = 246ms; resolution = 120,000), 200um (time = 118ms; resolution = 60,000).
- **Data analysis**: MaxQuant (v1.5.3.30) + Perseus (v1.5.6.0) + Origin Pro 2017

## Data collection

The data were collected from the PRIDE repository (accession ID: PXD008844). We downloaded the MaxQuant\_Peptides.txt and the MaxQuant\_ProteinGroups.txt files containing the combined identification and quantification results. The sample annotations were inferred from the names of columns holding the quantification data and the information in the article. The peptides data were converted to a `SingleCellExperiment::SingleCellExperiment` object. We split the protein table to separate the three types of quantification: protein intensity, label-free quantification (LFQ) and intensity based absolute quantification (iBAQ). Each table is converted to a `SingleCellExperiment::SingleCellExperiment` object along with the remaining protein annotations. The 4 objects are combined in a single `QFeatures::QFeatures` object and feature links are created based on the peptide leading razor protein ID and the protein ID.

## Source

The PSM data can be downloaded from the PRIDE repository PXD008844. FTP link <ftp://ftp.pride.ebi.ac.uk/pride/data/>

## References

Zhu, Ying, Maowei Dou, Paul D. Piehowski, Yiran Liang, Fangjun Wang, Rosalie K. Chu, William B. Chrisler, et al. 2018. “Spatially Resolved Proteome Mapping of Laser Capture Microdissected Tissue with Automated Sample Transfer to Nanodroplets.” *Molecular & Cellular Proteomics: MCP* 17 (9): 1864–74 ([link to article](#)).

## Examples

```
zhu2018MCP()
```

---

```
zhu2018NC_hela
```

```
Zhu et al. 2018 (Nat. Comm.): HeLa titration
```

---

## Description

Near single-cell proteomics data of HeLa samples containing different number of cells. There are three groups of cell concentrations: low (10-14 cells), medium (35-45 cells) and high (137-141 cells). The data also contain measures for blanks, HeLa lysates (50 cell equivalent) and 2 cancer cell line lysates (MCF7 and THP1, 50 cell equivalent).

## Usage

```
zhu2018NC_hela
```

## Format

A `QFeatures::QFeatures` object with 4 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `peptides`: quantitative information for 37,795 peptides from 21 samples
- `proteins_intensity`: protein intensities for 3,984 proteins from 21 samples
- `proteins_LFQ`: LFQ intensities for 3,984 proteins from 21 samples
- `proteins_iBAQ`: iBAQ values for 3,984 proteins from 21 samples

Sample annotation is stored in `colData(zhu2018NC_hela())`.

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the original article (see References).

- **Cell isolation:** HeLa cell concentration was adjusted by serial dilution and cell counting was performed manually using an inverted microscope.
- **Sample preparation** performed using the nanoPOTs device. Protein extraction using RapiGest (+ DTT) + alkylation (IAA) + Lys-C digestion + cleave RapiGest (formic acid).
- **Separation:** nanoACQUITY UPLC pump (60nL/min) with an Self-Pack PicoFrit 70cm x 30um LC columns.
- **Ionization:** ESI (1,900V).
- **Mass spectrometry:** Thermo Fisher Orbitrap Fusion Lumos Tribrid. MS1 settings: accumulation time = 246ms; resolution = 120,000; AGC = 1E6. MS/MS settings, depend on the sample size, excepted for the AGC = 1E5. Blank and approx. 10 cells (time = 502ms; resolution = 240,000), approx. 40 cells (time = 246ms; resolution = 120,000), approx. 140 cells (time = 118ms; resolution = 60,000).
- **Data analysis:** MaxQuant (v1.5.3.30) + Perseus + OriginLab 2017

## Data collection

The data were collected from the PRIDE repository (accession ID: PXD006847). We downloaded the CulturedCells\_peptides.txt and the CulturedCells\_proteinGroups.txt files containing the combined identification and quantification results. The sample annotations were inferred from the names of columns holding the quantification data and the information in the article. The peptides data were converted to a `SingleCellExperiment::SingleCellExperiment` object. We split the protein table to separate the three types of quantification: protein intensity, label-free quantification (LFQ) and intensity based absolute quantification (iBAQ). Each table is converted to a `SingleCellExperiment::SingleCellExperiment` object along with the remaining protein annotations. The 4 objects are combined in a single `QFeatures::QFeatures` object and feature links are created based on the peptide leading razor protein ID and the protein ID.

## Source

The PSM data can be downloaded from the PRIDE repository PXD006847. FTP link: <ftp://ftp.pride.ebi.ac.uk/pride/data/>

## References

Zhu, Ying, Paul D. Piehowski, Rui Zhao, Jing Chen, Yufeng Shen, Ronald J. Moore, Anil K. Shukla, et al. 2018. "Nanodroplet Processing Platform for Deep and Quantitative Proteome Profiling of 10-100 Mammalian Cells." Nature Communications 9 (1): 882 ([link to article](#)).

## See Also

The same experiment was conducted on HeLa lysates: [zhu2018NC\\_lysat](#).

## Examples

```
zhu2018NC_hela()
```

---

`zhu2018NC_islets`*Zhu et al. 2018 (Nat. Comm.): human pancreatic islets*

---

## Description

Near single-cell proteomics data human pancreas samples. The samples were collected from pancreatic tissue slices using laser dissection. The pancreata were obtained from organ donors through the JDRFNetwork for Pancreatic Organ Donors with Diabetes (nPOD) program. The sample come either from control patients (n=9) or from type 1 diabetes (T1D) patients (n=9).

## Usage

`zhu2018NC_islets`

## Format

A `QFeatures::QFeatures` object with 4 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `peptides`: quantitative information for 24,321 peptides from 18 islet samples
- `proteins_intensity`: quantitative information for 3,278 proteins from 18 islet samples
- `proteins_LFQ`: LFQ intensities for 3,278 proteins from 18 islet samples
- `proteins_iBAQ`: iBAQ values for 3,278 proteins from 18 islet samples

Sample annotation is stored in `colData(zhu2018NC_islets())`.

## Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation**: The islets were extracted from the pancreatic tissues using laser-capture microdissection.
- **Sample preparation** performed using the nanoPOTs device. Protein extraction using RapiGest (+ DTT) + alkylation (IAA) + Lys-C digestion + cleave RapiGest (formic acid)
- **Separation**: nanoACQUITY UPLC pump with an Self-Pack PicoFrit 70cm x 30um LC columns; 60nL/min)
- **Ionization**: ESI (1,900V)
- **Mass spectrometry**: Thermo Fisher Orbitrap Fusion Lumos Tribrid. MS1 settings: accumulation time = 246ms; resolution = 120,000; AGC = 1E6. MS/MS settings: accumulation time = 118ms; resolution = 60,000; AGC = 1E5.
- **Data analysis**: MaxQuant (v1.5.3.30) + Perseus + OriginLab 2017

## Data collection

The data were collected from the PRIDE repository (accession ID: PXD006847). We downloaded the `Islet_t1d_ct_peptides.txt` and the `Islet_t1d_ct_proteinGroups.txt` files containing the combined identification and quantification results. The sample types were inferred from the names of columns holding the quantification data. The peptides data were converted to a `SingleCellExperiment::SingleCellExperiment` object. We split the protein table to separate the three types

of quantification: protein intensity, label-free quantification (LFQ) and intensity based absolute quantification (iBAQ). Each table is converted to a `SingleCellExperiment::SingleCellExperiment` object along with the remaining protein annotations. The 4 objects are combined in a single `QFeatures::QFeatures` object and feature links are created based on the peptide leading razor protein ID and the protein ID.

### Source

The PSM data can be downloaded from the PRIDE repository PXD006847. The source link is: <ftp://ftp.pride.ebi.ac.uk/pride/data/archive/2018/01/PXD006847>

### References

Zhu, Ying, Paul D. Piehowski, Rui Zhao, Jing Chen, Yufeng Shen, Ronald J. Moore, Anil K. Shukla, et al. 2018. "Nanodroplet Processing Platform for Deep and Quantitative Proteome Profiling of 10-100 Mammalian Cells." *Nature Communications* 9 (1): 882 ([link to article](#)).

### Examples

```
zhu2018NC_islets()
```

---

zhu2018NC_lysates	<i>Zhu et al. 2018 (Nat. Comm.): HeLa lysates</i>
-------------------	---

---

### Description

Near single-cell proteomics data of HeLa lysates at different concentrations (10, 40 and 140 cell equivalent). Each concentration is acquired in triplicate.

### Usage

```
zhu2018NC_lysates
```

### Format

A `QFeatures::QFeatures` object with 4 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `peptides`: quantitative information for 14,921 peptides from 9 lysate samples
- `proteins_intensity`: quantitative information for 2,199 proteins from 9 lysate samples
- `proteins_LFQ`: LFQ intensities for 2,199 proteins from 9 lysate samples
- `proteins_iBAQ`: iBAQ values for 2,199 proteins from 9 lysate samples

Sample annotation is stored in `colData(zhu2018NC_lysates())`.

### Acquisition protocol

The data were acquired using the following setup. More information can be found in the original article (see References).

- **Cell isolation:** HeLas were collected from cell cultures.
- **Sample preparation** performed in bulk (5E5 cells/mL). Protein extraction using RapiGest (+ DTT) + dilution to target concentration + alkylation (IAA) + Lys-C digestion + trypsin digestion + cleave RapiGest (formic acid).
- **Separation:** nanoACQUITY UPLC pump (60nL/min) with an Self-Pack PicoFrit 70cm x 30um LC columns.
- **Ionization:** ESI (1,900V).
- **Mass spectrometry:** Thermo Fisher Orbitrap Fusion Lumos Tribrid. MS1 settings: accumulation time = 246ms; resolution = 120,000; AGC = 1E6. MS/MS settings, depend on the sample size, excepted for the AGC = 1E5. Blank and approx. 10 cells (time = 502ms; resolution = 240,000), approx. 40 cells (time = 246ms; resolution = 120,000), approx. 140 cells (time = 118ms; resolution = 60,000).
- **Data analysis:** MaxQuant (v1.5.3.30) + Perseus + OriginLab 2017.

### Data collection

The data were collected from the PRIDE repository (accession ID: PXD006847). We downloaded the `Vail_Prep_Vail_peptides.txt` and the `Vail_Prep_Vail_proteinGroups.txt` files containing the combined identification and quantification results. The sample annotations were inferred from the names of columns holding the quantification data and the information in the article. The peptides data were converted to a `SingleCellExperiment::SingleCellExperiment` object. We split the protein table to separate the three types of quantification: protein intensity, label-free quantification (LFQ) and intensity based absolute quantification (iBAQ). Each table is converted to a `SingleCellExperiment::SingleCellExperiment` object along with the remaining protein annotations. The 4 objects are combined in a single `QFeatures::QFeatures` object and feature links are created based on the peptide leading razor protein ID and the protein ID.

### Source

The PSM data can be downloaded from the PRIDE repository PXD006847. The source link is: <ftp://ftp.pride.ebi.ac.uk/pride/data/archive/2018/01/PXD006847>

### References

Zhu, Ying, Paul D. Piehowski, Rui Zhao, Jing Chen, Yufeng Shen, Ronald J. Moore, Anil K. Shukla, et al. 2018. "Nanodroplet Processing Platform for Deep and Quantitative Proteome Profiling of 10-100 Mammalian Cells." *Nature Communications* 9 (1): 882 ([link to article](#)).

### See Also

The same experiment was conducted directly on HeLa cells samples rather than lysates. The data is available in [zhu2018NC\\_hela](#).

### Examples

```
zhu2018NC_lysat()
```

---

zhu2019EL

*Zhu et al. 2019 (eLife): chicken utricle cells*

---

### Description

Single-cell proteomics data from chicken utricle acquired to study the hair-cell development. The cells are isolated from peeled utricular epithelium and separated into hair cells (FM1-43 high) and supporting cells (FM1-43 low). The sample contain either 1 cell (n = 28), 3 cells (n = 7), 5 cells (n = 8) or 20 cells (n = 14).

### Usage

zhu2019EL

### Format

A `QFeatures::QFeatures` object with 62 assays, each assay being a `SingleCellExperiment::SingleCellExperiment` object:

- `XYZw`: 60 assays containing PSM data. The sample are annotated as follows. X indicates the experiment, either 1 or 2. Y indicated the FM1-43 signal, either high (H) or low (L). Z indicates the number of cells (0, 1, 3, 5 or 20). w indicates the replicate, starting from a, it can go up to j.
- `peptides`: quantitative data for 3444 peptides in 60 samples (all runs are combined).
- `proteins_intensity`: protein intensities for 840 proteins from 24 samples
- `proteins_iBAQ`: iBAQ values for 840 proteins from 24 samples

Sample annotation is stored in `colData(zhu2019EL())`.

### Acquisition protocol

The data were acquired using the following setup. More information can be found in the source article (see References).

- **Cell isolation**: The cells were taken from the utricles of E15 chick embryos. Samples were stained with FM1-43FX and the cells were dissociated using enzymatic digestion. Cells were FACS sorted (BD Influx) and split based on their FM1-43 signal, while ensuring no debris, doublets or dead cells are retained.
- **Sample preparation** performed using the nanoPOTs device. Cell lysis and protein extraction and reduction are performed using dodecyl beta-D-maltoside + DTT + ammonium bicarbonate. Protein were then alkylated using IAA. Protein digestion is performed using Lys-C and trypsin. Finally samples acidification is performed using formic acid.
- **Separation**: Dionex UltiMate pump with an C18-Packed column (50cm x 30um; 60nL/min)
- **Ionization**: ESI (2,000V)
- **Mass spectrometry**: Orbitrap Fusion Lumos Tribrid. MS1 settings: accumulation time = 246ms; resolution = 120,000; AGC = 3E6. MS/MS settings: accumulation time = 502ms; resolution = 120,000; AGC = 2E5.
- **Data analysis**: Andromeda & MaxQuant (v1.5.3.30) and the search database is NCBI GRCg6a.

## Data collection

All data were collected from the PRIDE repository (accession ID: PXD014256).

The sample annotation information is provided in the `Zhu_2019_chick_single_cell_samples_CORRECTED.xlsx` file. This file was given during a personal discussion and is a corrected version of the annotation table available on the PRIDE repository.

The PSM data were found in the `evidence.txt` (in the `Experiment 1+ 2`) folder. The PSM data were filtered so that it contains only samples that are annotated. The data were then converted to a `QFeatures::QFeatures` object using the `scp::readSCP()` function.

The peptide data were found in the `peptides.txt` file. The column names holding the quantitative data were adapted to match the sample names in the `QFeatures::QFeatures` object. The data were then converted to a `SingleCellExperiment::SingleCellExperiment` object and then inserted in the `QFeatures::QFeatures` object. Links between the PSMs and the peptides were added

A similar procedure was applied to the protein data. The data were found in the `proteinGroups.txt` file. We split the protein table to separate the two types of quantification: summed intensity and intensity based absolute quantification (iBAQ). Both tables are converted to `SingleCellExperiment::SingleCellExperiment` objects and are added to the `QFeatures::QFeatures` object as well as the `AssayLink` between peptides and proteins.

## Source

The PSM data can be downloaded from the PRIDE repository PXD014256. The source link is: <ftp://ftp.pride.ebi.ac.uk/pride/data/archive/2019/11/PXD014256>

## References

Zhu, Ying, Mirko Scheibinger, Daniel Christian Ellwanger, Jocelyn F. Krey, Dongseok Choi, Ryan T. Kelly, Stefan Heller, and Peter G. Barr-Gillespie. 2019. "Single-Cell Proteomics Reveals Changes in Expression during Hair-Cell Development." *eLife* 8 (November). ([link to article](#)).

## Examples

`zhu2019EL()`

# Index

## \* datasets

brunner2022, [3](#)  
cong2020AC, [5](#)  
derks2022, [6](#)  
dou2019\_boosting, [8](#)  
dou2019\_lystates, [10](#)  
dou2019\_mouse, [12](#)  
gregoire2023\_mixCTRL, [14](#)  
guise2024, [16](#)  
hu2023\_K562, [18](#)  
khan2023, [20](#)  
krull2024, [22](#)  
leduc2022\_plexDIA, [24](#)  
leduc2022\_pSCoPE, [26](#)  
liang2020\_hela, [28](#)  
petrosius2023\_AstralAML, [29](#)  
petrosius2023\_mES, [31](#)  
schoof2021, [33](#)  
specht2019v2, [36](#)  
specht2019v3, [38](#)  
williams2020\_lfq, [40](#)  
williams2020\_tmt, [42](#)  
woo2022\_lung, [43](#)  
woo2022\_macrophage, [45](#)  
wu2026, [46](#)  
zhu2018MCP, [48](#)  
zhu2018NC\_hela, [49](#)  
zhu2018NC\_islets, [51](#)  
zhu2018NC\_lystates, [52](#)  
zhu2019EL, [54](#)

## \* dataset

ai2025a, [2](#)  
  
ai2025a, [2](#)  
  
brunner2022, [3](#)  
  
cong2020AC, [5](#)  
  
derks2022, [6](#)  
dou2019\_boosting, [8](#), [12](#), [14](#)  
dou2019\_lystates, [10](#), [10](#), [14](#)  
dou2019\_mouse, [10](#), [12](#), [12](#)  
  
gregoire2023\_mixCTRL, [14](#)

guise2024, [16](#)

hu2023 (hu2023\_K562), [18](#)  
hu2023\_K562, [18](#)  
hu2023\_oocyte (hu2023\_K562), [18](#)

khan2023, [20](#)  
krull2024, [22](#)

leduc2022, [23](#)  
leduc2022\_plexDIA, [24](#), [28](#)  
leduc2022\_pSCoPE, [25](#), [26](#)  
liang2020\_hela, [28](#)

petrosius2023\_AstralAML, [29](#)  
petrosius2023\_mES, [31](#)

QFeatures::addAssayLink, [9](#), [11](#), [13](#)  
QFeatures::AssayLinks, [41](#), [43](#), [44](#), [46](#)  
QFeatures::joinAssays, [9](#), [11](#), [13](#)  
QFeatures::QFeatures, [3–17](#), [19–46](#), [48–55](#)

schoof2021, [33](#)  
scp::readSCP(), [6](#), [9](#), [11](#), [13](#), [15](#), [17](#), [21](#), [27](#),  
[29](#), [31](#), [33](#), [34](#), [37](#), [39](#), [47](#), [55](#)

scpdata, [35](#)  
scpdata-package (scpdata), [35](#)  
SingleCellExperiment::SingleCellExperiment,  
[3–14](#), [16–55](#)

specht2019v2, [36](#)  
specht2019v3, [38](#)

williams2020\_lfq, [40](#)  
williams2020\_tmt, [42](#)  
woo2022\_lung, [43](#)  
woo2022\_macrophage, [45](#)  
wu2026, [46](#)

zhu2018MCP, [48](#)  
zhu2018NC\_hela, [49](#), [53](#)  
zhu2018NC\_islets, [51](#)  
zhu2018NC\_lystates, [50](#), [52](#)  
zhu2019EL, [54](#)