

# Package ‘tofsimsData’

May 19, 2026

**Type** Package

**Title** Import, process and analysis of ToF-SIMS imaging data

**Version** 1.41.0

**Date** 2014-10-23

**Author** Lorenz Gerber, Viet Mai Hoang

**Maintainer** Lorenz Gerber <genfys@gmail.com>

**Depends** R (>= 3.2.0)

**Description** This packages contains data to be used with the 'tofsims' package.

**License** GPL-3

**Suggests** knitr, rmarkdown, tools

**VignetteBuilder** knitr

**biocViews** ExperimentData, MassSpectrometry, ImagingMassSpectrometry,  
DataImport

**NeedsCompilation** no

**git\_url** <https://git.bioconductor.org/packages/tofsimsData>

**git\_branch** devel

**git\_last\_commit** 0a2ebb4

**git\_last\_commit\_date** 2026-04-28

**Repository** Bioconductor 3.24

**Date/Publication** 2026-05-19

## Contents

tofsimsData-package . . . . .	2
testImage . . . . .	2
testSpectra . . . . .	3
<b>Index</b>	<b>4</b>

---

tofsimsData-package    *tofsimsData*

---

### Description

ToF-SIMS Toolbox

### Details

Package:    tofsimsData  
Type:        Package  
Version:    0.99  
Date:        23-10-2014  
License:    GPL (>=2)  
LazyLoad:   yes

Toolbox for Time-of-Flight Secondary Ion Mass-Spectrometry (ToF-SIMS) data processing and analysis. The package facilitates importing of raw data files, loading preprocessed data and a range of multivariate analysis methods that are most commonly applied in imaging (ToF-SIMS) mass spectrometry.

### Author(s)

Lorenz Gerber <lorenz.gerber@slu.se>

---

testImage                    *Example ToF-SIMS data*

---

### Description

A dataset containing a MassImage recorded on a Ulvac-Phi TRIFT-II ToF-SIMS. The .RAW data file was imported using `tofsimsImage<-MassImage('ulvacrawpeaks','filename',PeakList=tofsimsSpectra)`. The sample is a freeze-dried transversal poplar wood section of 100 micrometer thickness.

### Usage

```
data(tofsimsData)
```

### Format

A MassImage object

### Value

MassImage object

---

`testSpectra`*Example ToF-SIMS data*

---

**Description**

A dataset containing a `MassSpectra` recorded on a Ulvac-Phi TRIFT-II ToF-SIMS. The `.RAW` data file was imported using `tofsimsSpectra<-MassSpectra('ulvacraw', 'filename')`. The sample is a freeze-dried transversal poplar wood section of 100 micrometer thickness.

**Usage**

```
data(tofsimsData)
```

**Format**

A `MassSpectra` object

**Value**

`MassSpectra` object

# Index

\* **dataset**

testImage, [2](#)  
testSpectra, [3](#)

\* **package**

tofsimsData-package, [2](#)

testImage, [2](#)

testSpectra, [3](#)

tofsimsData-package, [2](#)